



SKU:

## **HERO MEDICAL OUTDOOR LITE SEAT WALKER (ALUMINIUM)**

**\$499.99**

Exclusive to ILS, the Hero Medical Outdoor Lite Seat Walker (Aluminium) blends functionality with comfort for outdoor independence. Available in red or black and three sizes, it offers ergonomic design and secure mobility in a lightweight frame.



Scan QR code to  
view product online

## VARIATIONS





Image	SKU	Price	Description	Height	Colours
	S301678	\$499.99		Tall	Black
	S301675	\$499.99		Tall	Red
	S301674	\$499.99		Standard	Red
	S301677	\$499.99		Standard	Black

Image	SKU	Price	Description	Height	Colours
	S301676	\$499.99		Low	Red
	S301679	\$499.99		Low	Black

## PRODUCT DESCRIPTION

The Hero Medical Outdoor Lite Seat Walker (Aluminium) is designed to promote stable, comfortable movement outdoors. With a curved lightweight aluminium frame (6.4 kg – 6.6 kg), it's easy to manoeuvre while offering reliable support.

Key features include large 9.5" front castors for stable travel, light-touch handbrakes with a parking function, and a padded, height-adjustable backrest with four clip settings. The spill-resistant nylon seat and waterproof storage bag provides comfort and practicality.

This model folds quickly for storage and transport and is delivered pre-assembled for immediate use.

### Who it's for

- For clients seeking a stylish, lightweight walker for confident, daily independent mobility.
- For therapists prescribing a reliable aluminium walking aid that supports posture, comfort and safe mobility.

### Key features

- **Lightweight aluminium frame (6.4 kg – 6.6 kg)** – stable yet easy to handle
- **Large 9.5" front castors** – support smooth outdoor manoeuvrability
- **Height-adjustable rounded backrest** – four clip settings for individualised comfort
- **Ergonomic height-adjustable handles** – with quick adjustments for custom fit
- **Light-touch handbrakes** – with integrated parking brake
- **Foam-padded, spill-resistant seat** – supports rest and comfort
- **Waterproof nylon storage bag** – removable and easy to access
- **Simple folding mechanism** – enables compact storage and portability
- **Delivered pre-assembled** – ready for use after unfolding

## ADDITIONAL INFORMATION

<b>Brand</b>	Hero Medical
<b>Max safe working load (SWL)</b>	136 kg (299 lbs)
<b>Overall height</b>	Low - 82.5 - 95 cm (32.4 - 37.4 inches) Standard - 84.5 - 97 cm (33.2 - 38.1 inches) Tall - 93.5 - 106 cm (36.8 - 41.7 inches)
<b>Seat height</b>	Low - 50.5 cm (19.8 inches) Standard - 52.5 cm (20.6 inches) Tall - 60.5 cm (23.8 inches)
<b>Seat width</b>	Low - 46 cm (18.1 inches) Standard - 46 cm (18.1 inches) Tall - 46 cm (18.1 inches)
<b>Outer width (including back castors)</b>	Low - 60 cm (23.6 inches) Standard - 60 cm (23.6 inches) Tall - 60 cm (23.6 inches)
<b>Inner width (between back frames)</b>	Low - 46.5 cm (18.3 inches) Standard - 46.5 cm (18.3 inches) Tall - 45.5 cm (18.3 inches)
<b>Product frame weight</b>	Low - 6.4 kg (14.1 lbs) Standard - 6.5 kg (14.3 lbs) Tall - 6.6 kg (14.5 lbs)
<b>Gross weight</b>	Low - 7.8 kg (17.1 lbs) Standard - 7.9 kg (17.4 lbs) Tall - 8 kg (17.6 lbs)
<b>Frame material</b>	Alumimium
<b>Colour options</b>	Black and red
<b>Warranty period</b>	12 months for structural frame