




**SKU:**

## **HERO MEDICAL OUTDOOR CARBON FIBRE LITE WALKER**

**\$899.99**

Exclusive to ILS, the Hero Medical Outdoor Carbon Fibre Lite Walker has been designed to maximise your mobility and independence. The outdoor lite walker has many premium luxury features to ensure you experience effortless mobility outdoors including an ultra-lightweight carbon fibre frame of 5.7 kg, large front wheels, sturdy braking and safety systems and easy-to-operate ergonomic functions.



Scan QR code to  
view product online

**VARIATIONS**

| Image  | SKU     | Price    | Description | Height   |
|--|---------|----------|-------------|----------|
|   | S21707  | \$899.99 |             | Standard |
|  | S301673 | \$899.99 |             | Tall     |

## PRODUCT DESCRIPTION

### **Hero Medical Outdoor Lite Walker - Carbon Fibre Walking Frame**

The Hero Medical Outdoor Carbon Fibre Lite Walker offers premium luxury features to ensure you experience effortless mobility outdoors. As a heavy-duty quad seat walker, it has a safe working limit of 136 kg. The Hero Medical Outdoor Lite Walker is the perfect blend of lightness and superior durability. Boasting a unique curved design, the Hero Medical Outdoor Lite Walker is super stylish and available in two sizes (standard and tall).

#### **Key features:**

- Ultra lightweight - carbon fibre frame is only 5.7 kg
- Large 9.5" front castors for easy maneuvering and stability
- Padded easy-fold rounded adjustable backrest for comfortable seating. There are four levels of adjustment which can be easily made by pressing a clip.
- Adjustable handle height - 6 levels of adjustment
- Light touch hand brake system with an integrated parking brake function
- Super easy to fold and store
- Removable stylish tote storage bag that straps easily on the seat walker frame
- Comfortable padded seat with integrated layer of foam
- Ergonomic and adjustable handle bars with height adjustment clips
- Comes assembled in the box and just needs to be unfolded

#### **Free bonus accessory:**

- Walking stick/cane holder which easily clips on the seat walker frame

## ADDITIONAL INFORMATION

|                             |  |
|-----------------------------|--|
| <b>Brand</b>                | Hero Medical   |
| <b>Max SWL</b>              | 136 kg (299 lbs)   |
| <b>Overall height</b>       | Standard – 84.5 – 97cm (33.2-38.1 inches)<br>Tall 93.5 – 106 cm (36.8 – 41.7 inches) |
| <b>Seat height</b>          | Standard 52.5 cm (20.6 inches)<br>Tall 60.5 cm (23.8 inches)                         |
| <b>Seating width</b>        | Standard 45 cm (17.7 inches)<br>Tall 45 cm (17.7 inches)                             |
| <b>Outer width</b>          | Standard 58 cm (22.8 inches)<br>Tall 58 cm (22.8 inches)                             |
| <b>Inner width</b>          | Standard 45.5 cm (17.9 inches)<br>Tall 45.5 cm (17.9 inches)                         |
| <b>Product frame weight</b> | Standard 5.7 kg (12.5 lbs)<br>Tall 5.9 kg (13 lbs)                                   |
| <b>Gross weight</b>         | Standard 7.1 kg (15.6 lbs)<br>Tall 7.3 kg (16 lbs)                                   |
| <b>Frame material</b>       | Carbon fibre   |
| <b>Colour</b>               | Black  |
| <b>Warranty period</b>      | 12 months for structural frame   |