



SKU:

HERO MEDICAL OUTDOOR LITE SEAT WALKER (CARBON FIBRE)

\$899.99

Exclusive to ILS, the Hero Medical Outdoor Lite Seat Walker (Carbon Fibre), provides ergonomic walking support with quality materials and a curved frame design. Available in Standard and Tall sizes, it's ideal for those who need supreme performance and durability.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Height
	S21707	\$899.99		Standard
	S301673	\$899.99		Tall

PRODUCT DESCRIPTION

The Hero Medical Outdoor Lite Seat Walker (Carbon Fibre) features a curved carbon fibre frame weighing just 5.7 – 5.9 kg, providing effortless manoeuvrability without compromising strength. Designed for daily indoor and outdoor use, it combines stability, comfort and style.

The large 9.5” front castors handle varied terrain with ease.

The walker features a comfortable padded seat and rounded backrest, with four clip-adjustable heights, to provide support during rest.

The handle height adjustment provides ergonomic walking posture.

This walker includes a stylish removable tote bag and folds compactly for easy transport.

A free cane holder accessory also clips conveniently onto the frame.

Who it's for

- For clients seeking an ultra-lightweight, high performance seat walker for indoor and outdoor independence.
- For therapists prescribing an ultra-lightweight, high-performance seat walker that combines advanced functionality, long-term durability, and contemporary design.

Key features

- **Ultra-light carbon fibre frame (5.7 – 5.9 kg)** – supports easy lifting and transport
- **Large 9.5” front castors** – smooth travel across indoor and outdoor surfaces
- **Height-adjustable rounded backrest** – four clip settings for comfort
- **Ergonomic height-adjustable handles** – with quick adjustments for custom fit
- **Light-touch handbrakes** – with integrated parking brake
- **Foam-padded seat** – offers supportive, durable rest surface
- **Removable tote-style storage bag** – stylish and practical for personal items
- **Quick-folding frame** – enables compact transport and storage
- **Delivered pre-assembled** – ready for immediate use
- **Free cane holder accessory** – clips onto frame for added convenience

ADDITIONAL INFORMATION

Brand	Hero Medical
Max SWL	136 kg (299 lbs)
Overall height	Standard – 84.5 – 97cm (33.2-38.1 inches) Tall 93.5 – 106 cm (36.8 – 41.7 inches)
Seat height	Standard 52.5 cm (20.6 inches) Tall 60.5 cm (23.8 inches)
Seating width	Standard 45 cm (17.7 inches) Tall 45 cm (17.7 inches)
Outer width	Standard 58 cm (22.8 inches) Tall 58 cm (22.8 inches)
Inner width	Standard 45.5 cm (17.9 inches) Tall 45.5 cm (17.9 inches)
Product frame weight	Standard 5.7 kg (12.5 lbs) Tall 5.9 kg (13 lbs)
Gross weight	Standard 7.1 kg (15.6 lbs) Tall 7.3 kg (16 lbs)
Frame material	Carbon fibre
Colour	Black
Warranty period	12 months for structural frame