



SKU:


HERO MEDICAL SUPER BUDGET WHEELCHAIR

\$349.99

The Hero Medical Super Budget Wheelchair, available only at ILS, is a robust, entry-level wheelchair that offers exceptional value, and is available in both self-propelled and transit options.

Crafted with a sturdy steel frame and durable, easy to clean vinyl upholstery, this wheelchair features an 18" seat and effortless folding capability for convenient transport and storage.

Both models come equipped with essential safety features, including park brakes, seatbelt, and height-adjustable swing-away legrests, ensuring user comfort and security.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Type
	S310064	\$349.99		Self-propelled
	S310065	\$349.99		Transit

PRODUCT DESCRIPTION

The Super Budget Wheelchair offers robustness and exceptional value, combining durability with affordability. Both the self-propelled and transit options are practical, easy to maintain and are suitable for use within both community and facility settings. Both models come with a seat width of 18" and a safe work load of 115 kg.

Key Features:

- **Sturdy:** Built with a robust and durable steel frame.
- **Easy to manage and clean:** Features vinyl material on the seat and backrest for effortless maintenance, ideal for facility settings. Folds easily for transportation.
- **Comfortable:** Padded seat, backrest, and armrests and height adjustable footplates.
- **Great Value:** Offers affordability without compromising on features and benefits.
- **Safety First:** Equipped with essential safety features, including a park brake, seatbelt, anti-tippers, and swing-away leg rests. Compliant with AS/NZS 3695.1:2011

ADDITIONAL INFORMATION

Brand	Hero Medical
Weight capacity	115kg (253 lbs)
Length	106cm (42")
Height	90cm (36")
Seat width	45cm (18")
Seat length	40cm (16")
Seat height	49cm (19")
Unit weight	18kg (39 lbs)
Type	Self-propelled, Transit