


SKU:

HIGHBACK BC1 SUPER KINGSIZE BARIATRIC CHAIR

\$1,299.99

The Highback BC1 Super Kingsize Bariatric Chair is an ergonomically designed single-piece chair with adjustable leg height.



Scan QR code to
view product online



VARIATIONS

Image	SKU	Price	Description	Colour
	S10308	\$1,299.99		Latte
	S10306	\$1,299.99		Blue

PRODUCT DESCRIPTION

The Highback BC1 Super Kingsize Bariatric Chair is designed as an ergonomically single-piece shell with a steel frame that is fully welded and leg height that can be adjusted with foam that is fire retardant.

ADDITIONAL INFORMATION

Overall Height	99- 114 cm
Depth	74 cm
Safe Working Load	300 Kg
Seat width	74cm
Seat Height	36cm – 51cm
Back Height	72cm