



**SKU:**

## **DAY CHAIR BARIATRIC HIGHBACK SUPER KING BC2 300 KG**

**\$1,499.99**

The BC2 Super King size day chair is extremely durable and comfortable, featuring an ergonomic additional segmented backrest and seat pad to provide added comfort and functionality for its users. Perfectly suited for aged care, retirement, healthcare, hospital and homecare settings.



Scan QR code to  
view product online

## VARIATIONS

Image	SKU	Price	Description	Select Colour
	S9941	\$1,499.99		Blue Vinyl
	S10309	\$1,499.99		Latte Vinyl

## **PRODUCT DESCRIPTION**

The BC2 Super King Day chair consists of a fully welded powder coated tubular steel frame with height adjustable legs. These chairs come with ergonomic additional segmented backrest and seat pad for added comfort and ease of cleaning. High weight capacity of up to 300 kg! Vinyl upholstery available in blue and latte colours.

## **ADDITIONAL INFORMATION**

<b>Safe working load</b>	300 Kg (661 lbs)
<b>Upholstery</b>	Vinyl
<b>Warranty</b>	10 years for structural component and 5 years for steel frame