



SKU:

BC2 SUPER KINGSIZE BARIATRIC HIGH BACK DAY CHAIR

\$1,499.99

The BC2 Super Kingsize Bariatric Day Chair is a heavy-duty seating solution offering maximum comfort and strength for bariatric clients. With an extra-wide seat, segmented backrest and height-adjustable frame, it's suited to high-needs environments across aged care, clinical and home settings.



Scan QR code to
view product online

VARIATIONS

SKU	Price	Description	Colours
S9941	\$1,499.99	The BC2 Super Kingsize Bariatric Day Chair is a heavy-duty seating solution offering maximum comfort and strength for bariatric clients. With an extra-wide seat, segmented backrest and height-adjustable frame, it's suited to high-needs environments across aged care, clinical and home settings.	Blue Vinyl
S10309	\$1,499.99	The BC2 Super Kingsize Bariatric Day Chair is a heavy-duty seating solution offering maximum comfort and strength for bariatric clients. With an extra-wide seat, segmented backrest and height-adjustable frame, it's suited to high-needs environments across aged care, clinical and home settings.	Latte Vinyl

PRODUCT DESCRIPTION

The BC2 Super Kingsize Bariatric Day Chair delivers enhanced seating support for bariatric clients who require greater width, strength and seated comfort.

It features a fully welded, powder-coated tubular steel frame with a weight capacity of up to 300 kg, making it ideal for long-term use in high-demand care environments. The Super King model offers an extra-wide 74 cm seat and segmented backrest design to promote pressure relief and postural support.

The seat height can be adjusted between 41 cm - 56 cm to suit individual needs and support safer, more ergonomic transfers. The padded seat and back cushions are upholstered in easy-to-clean vinyl and available in blue or latte colour options.

Designed for both comfort and durability, it's well suited to aged care, hospital, rehabilitation and home use.

Who it's for

- For clients who require extra-wide, high-capacity seating with ergonomic comfort and support.
- For therapists seeking an adjustable, bariatric day chair for their client.

Key features

- **Extra-wide 74 cm seat** – supports larger clients requiring more space and comfort
- **300 kg weight capacity** – heavy-duty construction for high-care bariatric needs
- **Fully welded powder-coated steel frame** – ensures long-term strength and stability
- **Segmented padded backrest and seat** – improves pressure relief and comfort
- **Height-adjustable legs** – allows tailored seat height for safe and supported transfers
- **Vinyl upholstery** – available in blue or latte and easy to wipe clean

ADDITIONAL INFORMATION

Safe Working Load (SWL)	300 kg
Overall Width	98 cm
Overall Depth	74 cm
Overall Height	99 cm – 114 cm
Seat Width	74 cm
Seat Depth	49 cm
Seat to Floor Height	36 cm – 51 cm
Backrest Height	N/A
Brand	Atama Furniture

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

No FAQ Found