




**SKU:** S18675

## **INVACARE BIRDIE 180KG PATIENT LIFT HOISTS**

**From: \$4,499.99**

The Birdie Evo Hoist has been designed to offer a high feature-set of options that have been blended with new technologies to provide maximum comfort and security when lifting or transferring a patient to or from a bed, chair or even the floor. Users can be easily rotated 360°, making positioning much easier for carers, and reducing anxiety for clients. The leg width is electrically or mechanically adjustable, with a minimum width of 64cm and a maximum spread up to 102cm. The Birdie EVO is easily folded without needing the use of tools, thereby taking up minimal space & making storage & transport easier.



Scan QR code to  
view product online

## **PRODUCT DESCRIPTION**

- Maximised space: the high boom ensures there is space for the client, even when the hoist is in its highest position and the optimised space in front of the actuator reduces the risk of injury to the clients' knees
- Easy rotation: the user can be easily rotated 360° while being lifted, promoting easier transfers and positioning between bed, chair or wheelchair
- Improved ergonomic hook design for easier sling attachment. This design allows the loop on the sling to be attached with just one hand. A wider hook makes positioning loops easier, whilst the hook's nose has been designed to prevent slings sliding out unintentionally, reducing the risk of injury
- Optional 4-point spreader bar – accessory (S300789)
- Smartlock: a new quick release mechanism makes it simple and safe to exchange the spreader bar without tools. Enhanced safety is provided via a swallow tail connection that is further secured by a locking pin to ensure both parts of the system remain in place
- Slow'R: a new dampener has been designed and integrated into the boom to reduce the rocking movement of the client, giving them a more comfortable experience and providing carers with more control during lifts and transfers
- Tool-less folding for compact storage
- Electric leg adjustment
- Includes 2-point spreader bar
- Detachable battery pack
- Ergonomic push bar makes handling and manoeuvring even easier.

## ADDITIONAL INFORMATION

<b>Total Width Closed</b>	680 mm
<b>Total Length</b>	1150 mm
<b>Total Internal Width (open / closed)</b>	1040 mm / 580 mm
<b>Lifting Range (from attach point)</b>	470 mm – 1740 mm
<b>Turning Radius</b>	1400 mm
<b>Base Height (top of legs)</b>	115 mm
<b>Total weight (incl. leg spreader)</b>	44 kg
<b>Max SWL</b>	180 kg
<b>Brand</b>	Invacare
<b>Available on Queensland Medical Aids Subsidy Scheme</b>	Yes