



HOIST JASMINE PATIENT TRANSFER

~~\$4,195.99~~ **\$3,695.99**

The Jasmine Mobile Lift is a flexible solution that offers high comfort for all clients. The lift covers all needs, providing safe patient handling with up to 200kg weight capacity. With its expansive, wide opening base and high lifting mechanism, the Jasmine makes life easy for caregivers operating this Invacare lift system. The control box with diagnostic system is an intelligent service monitoring system that enables a safe operating environment for the caregiver, whilst allowing maximum comfort for the client during each transfer.



SKU: S8447

Categories: [Hoists](#), [Patient Transfer Equipment](#)

GALLERY IMAGES



PRODUCT DESCRIPTION

Power base legs adjust easily with push-button hand pendant. Total pinch-point protection covers all moving parts for added safety. Height range allows lifting from floor and high surfaces. Safety-tested for highest quality standards. Powered lift mechanism allows transfer of up to 500lbs with the push of the button.

ADDITIONAL INFORMATION

Brand	Invacare
Manufacturer's part number	1156008
Max SWL	200kg
Total Width Closed	690mm
Total Width Open	1050mm
Total Length	1340mm

Weight	43kg
Lifting Range	54cm - 192cm
Available on Queensland Medical Aids Subsidy Scheme	Yes