




SKU:

I-CARE COMPANION BED - RANGE





~~\$1,699.99~~ \$1,699.99

The I-Care Companion Base offers a stationary sleeping surface for partners or caregivers, compatible with I-Care adjustable homecare beds. It provides a supportive, non-clinical solution for shared spaces without compromising care.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Select Size	Select Colour
	S11280	\$1,499.99		Long Single	Black
	S13201	\$1,499.99		Long Single	Grey
	S11341	\$1,699.99		King Single	Black
	S13942	\$1,699.99		King Single	Grey

PRODUCT DESCRIPTION

The I-Care Companion Base supports partners, parents or caregivers who wish to sleep next to a loved one without interfering with their clinical care. It's a stationary bed designed for use with the I-Care homecare bed range.

The companion base has no adjustable functions, ensuring stability and simplicity. It features a timber slat sleeping surface and sits at a fixed height of 34.5 cm from the floor to the mattress platform. The lightweight frame and underbed clearance allow easy repositioning and cleaning access.

Finished in durable commercial-grade fabric, it matches seamlessly with I-Care bed models and is available in stone grey or onyx black upholstery.

Who it's for

- For partners, parents or caregivers wanting to sleep beside someone in an I-Care homecare bed.
- For therapists wanting to support family-centred care with a fixed-height companion base.

Key Features

- **Stationary bed frame** – supports shared sleep without clinical adjustment features
- **Fixed height platform** – sits at 34.5 cm to align with I-Care bed ranges
- **Timber slat base** – to ensure support and breathability
- **Lightweight and movable** – easy to reposition as needed
- **Underbed clearance** – allows cleaning access and cable routing
- **Commercial-grade fabric** – matches I-Care adjustable beds

ADDITIONAL INFORMATION

Brand	I-Care Medical
Select Size	Long Single, King Single
Select Colour	Black, Grey
Bed sizing	Long Single – 2030 x 900 mm King Single – 2030 x 1070 mm