



SKU:

## I-CARE IC25 MATTRESS





~~\$1,899.99~~ **\$1,899.99**

The I-Care IC25 visco elastic mattress is made of high density, open-cell material that responds to temperature, and provides pressure reduction.



Scan QR code to  
view product online

**VARIATIONS**

Image	SKU	Price	Description	Size
	S10585	\$899.99		Long Single
	S10584	\$999.99		King Single
	S305285	\$1,799.99		Long Double
	S10583	\$1,899.99		Queen

## PRODUCT DESCRIPTION

### Pressure Redistribution

ActiveX™ Material is always adapting to allow blood to circulate freely and to minimise pressure on areas of risk. It responds to pressure and temperature.

### No Springs

Icare mattresses don't have any springs which makes them suitable for any profiling or static mattress platforms. We have found springs cause problems over time particularly on medical beds. And of course springs don't have a good pressure care track record.

### Key features

- 22.5 cm Mattress Depth
- SWL Minimum 30 kg
- SWL Maximum 200 kg
- For use with any I-Care Bed
- The mattresses are also available in two more variations: medium support [IC20](#) and supportive firm [IC15](#)

## ADDITIONAL INFORMATION

<b>Brand</b>	I-Care Medical
<b>Softness</b>	Soft Plush
<b>Size</b>	Double, Long Double, King Single, Long Single, Queen
<b>Warranty</b>	10 Years
<b>Max SWL</b>	200 kg
<b>Long single mattress size</b>	203 x 90 cm
<b>King single mattress size</b>	203 x 107 cm
<b>Double mattress size</b>	190 x 135 cm