



SKU:

## I-CARE IC333 ULTRA- LO HOSPITAL BED

~~\$7,499.99~~ \$7,499.99

The IC333 Bed is unmatched for function and durability, making the hospital look out of your room. This four-function bed features dual lift motors totalling 12,000 newtons lifting power, unlimited under-bed clearance, and eight heavy-duty casters.

\*Does not come with mattress.



Scan QR code to  
view product online

**VARIATIONS**

<b>Image</b>	<b>SKU</b>	<b>Price</b>	<b>Description</b>	<b>Bed size</b>
	S10589	\$4,299.99		Long Single
	S9907	\$4,699.99		King Single
	S15617	\$6,999.99		Double
	S10590	\$7,499.99		Queen

## **PRODUCT DESCRIPTION**

Key features:

- Strong Breathable Mattress Platform
- Available With Headboards
- External Rechargeable Battery Backup
- Lockable Heavy Duty Castors
- End And Side Mattress Retainer Bars
- Back Lit Wipeable Hand Piece
- No Gap Between Mattresses On King Beds (2x Long Singles)
- Unlimited Under Bed Clearance
- Ultra Low Design Reducing Risk Of Falls

\*Does not come with mattress.

## **ADDITIONAL INFORMATION**

<b>Brand</b>	I-Care
<b>Max SWL</b>	Long Single: 200kg / 441 lbs SWL King Single: 200kg / 441 lbs SWL Double: 300kg / 661 lbs SWL Queen: 300kg / 661 lbs SWL
<b>Top of bed max height</b>	66 cm (25.9 inches)
<b>Top of bed min height</b>	22 cm (8.6 inches)
<b>Warranty</b>	Steel Frame of Bed - 10 Year Limited Warranty Actuators - 5 Year Limited Warranty Hand Piece - 2 Year Limited Warranty Control Box - 5 Year Limited Warranty Upholstered Surround and Fabric - 3 Year Limited Warranty Castor wheels - 5 Year Limited Warranty
<b>Long single sizing details</b>	Mattress Surface Size: 2030 x 900 Base Size: 2150 x 960
<b>King single sizing details</b>	Mattress Surface Size: 2030 x 1070 Base Size: 2150 x 1120
<b>Long double sizing details</b>	Mattress Surface Size: 2030 x 1350 Base Size: 2150 x 1400 One hand control
<b>Queen sizing details</b>	Mattress Surface Size: 2030 x 1500 Base Size: 2150 x 1580 One hand control