



SKU:

I-CARE IC333 ULTRA- LO HOSPITAL BED

~~\$9,999.99~~ \$9,999.99

The IC333 range of beds are an ideal functional and durable homecare bed solution. This four-function bed features dual lift motors totalling 12,000 newtons lifting power, unlimited under bed clearance and eight heavy duty castors.



Scan QR code to
view product online

VARIATIONS





Image	SKU	Price	Description	Colour	Size
	S302366	\$9,999.99		Grey	Split Queen
	S302365	\$9,999.99		Black	Split Queen
	S11578	\$4,299.99		Grey	Long Single
	S10589	\$4,299.99		Black	Long Single





Image	SKU	Price	Description	Colour	Size
	S11570	\$4,699.99		Grey	King Single
	S9907	\$4,699.99		Black	King Single
	S10785	\$6,999.99		Grey	Double
	S15617	\$6,999.99		Black	Double

Image	SKU	Price	Description	Colour	Size
	S16616	\$7,499.99		Grey	Queen
	S10590	\$7,499.99		Black	Queen

PRODUCT DESCRIPTION

- Head Lift/Backrest Raise Function
- Hi-Lo Height Adjustable Function
- Knee Break Function
- Trendelenburg and Reverse Trendelenburg Function
- Backlit Wipeable Hand Piece
- 200kg SWL on Long Single and King Single
- 300kg SWL on Double and Queen
- Lockable Castors
- Hard Wearing Commercial Grade Fabric
- Head Lift Angle 62°
- Knee Break Angle 42°
- Height 22cm-66cm (*Measured from floor to top of mattress platform*)
- Strong Breathable Mattress Platform
- End and Corner Mattress Retainer Bars
- Please note: Mattress and Accessories are Additional

ADDITIONAL INFORMATION

Brand	I-Care
Max SWL	Long Single: 200kg / 441 lbs SWL King Single: 200kg / 441 lbs SWL Double: 300kg / 661 lbs SWL Queen: 300kg / 661 lbs SWL
Top of bed max height	66 cm (25.9 inches)
Top of bed min height	22 cm (8.6 inches)
Colour	Black, Grey
Warranty	Steel Frame of Bed – 10 Year Limited Warranty Actuators – 5 Year Limited Warranty Hand Piece – 2 Year Limited Warranty Control Box – 5 Year Limited Warranty Upholstered Surround and Fabric – 3 Year Limited Warranty Castor wheels – 5 Year Limited Warranty
Long single sizing details	Mattress Surface Size: 2030 x 900 Base Size: 2150 x 960
King single sizing details	Mattress Surface Size: 2030 x 1070 Base Size: 2150 x 1120
Long double sizing details	Mattress Surface Size: 2030 x 1350 Base Size: 2150 x 1400 One hand control
Queen sizing details	Mattress Surface Size: 2030 x 1500 Base Size: 2150 x 1580 One hand control