



SKU:

## I-CARE IC333 ULTRA- LO HOSPITAL BED





~~\$9,999.99~~ \$9,999.99





The IC333 range of beds are an ideal functional and durable homecare bed solution. This four-function bed features dual lift motors totalling 12,000 newtons lifting power, unlimited under bed clearance and eight heavy duty castors.



Scan QR code to  
view product online

**VARIATIONS**

Image	SKU	Price	Description	Colour	Size
	S302366	\$9,999.99		Grey	Split Queen
	S302365	\$9,999.99		Black	Split Queen
	S11578	\$4,299.99		Grey	Long Single
	S10589	\$4,299.99		Black	Long Single

<b>Image</b>	<b>SKU</b>	<b>Price</b>	<b>Description</b>	<b>Colour</b>	<b>Size</b>
	S11570	\$4,699.99		Grey	King Single
	S9907	\$4,699.99		Black	King Single
	S10785	\$6,999.99		Grey	Double
	S15617	\$6,999.99		Black	Double

<b>Image</b>	<b>SKU</b>	<b>Price</b>	<b>Description</b>	<b>Colour</b>	<b>Size</b>
	S16616	\$7,499.99		Grey	Queen
	S10590	\$7,499.99		Black	Queen

## **PRODUCT DESCRIPTION**

- Head Lift/Backrest Raise Function
- Hi-Lo Height Adjustable Function
- Knee Break Function
- Trendelenburg and Reverse Trendelenburg Function
- Backlit Wipeable Hand Piece
- 200kg SWL on Long Single and King Single
- 300kg SWL on Double and Queen
- Lockable Castors
- Hard Wearing Commercial Grade Fabric
- Head Lift Angle 62°
- Knee Break Angle 42°
- Height 22cm-66cm (*Measured from floor to top of mattress platform*)
- Strong Breathable Mattress Platform
- End and Corner Mattress Retainer Bars
- Please note: Mattress and Accessories are Additional

## **ADDITIONAL INFORMATION**

<b>Brand</b>	I-Care
<b>Max SWL</b>	Long Single: 200kg / 441 lbs SWL King Single: 200kg / 441 lbs SWL Double: 300kg / 661 lbs SWL Queen: 300kg / 661 lbs SWL
<b>Top of bed max height</b>	66 cm (25.9 inches)
<b>Top of bed min height</b>	22 cm (8.6 inches)
<b>Colour</b>	Black, Grey
<b>Warranty</b>	Steel Frame of Bed – 10 Year Limited Warranty Actuators – 5 Year Limited Warranty Hand Piece – 2 Year Limited Warranty Control Box – 5 Year Limited Warranty Upholstered Surround and Fabric – 3 Year Limited Warranty Castor wheels – 5 Year Limited Warranty
<b>Long single sizing details</b>	Mattress Surface Size: 2030 x 900 Base Size: 2150 x 960
<b>King single sizing details</b>	Mattress Surface Size: 2030 x 1070 Base Size: 2150 x 1120
<b>Long double sizing details</b>	Mattress Surface Size: 2030 x 1350 Base Size: 2150 x 1400 One hand control
<b>Queen sizing details</b>	Mattress Surface Size: 2030 x 1500 Base Size: 2150 x 1580 One hand control