



**SKU:**

## **ICARE IC333 HOMECARE BED - RANGE**





**~~\$9,999.99~~ \$9,999.99**





The Icare IC333 Homecare Bed combines clinical-grade functionality with a home-friendly design. Its ultra-low height, full adjustability and durable build make it a reliable option for long-term care at home.



Scan QR code to  
view product online

## VARIATIONS

| Image   | SKU     | Price      | Description | Colour | Sizes       |
|---|---------|------------|-------------|--------|-------------|
|     | S302366 | \$9,999.99 |             | Grey   | Split Queen |
|    | S302365 | \$9,999.99 |             | Black  | Split Queen |
|  | S11578  | \$4,299.99 |             | Grey   | Long Single |
|  | S10589  | \$4,299.99 |             | Black  | Long Single |

| Image  | SKU    | Price      | Description | Colour | Sizes       |
|--|--------|------------|-------------|--------|-------------|
|   | S11570 | \$4,699.99 |             | Grey   | King Single |
|   | S9907  | \$4,699.99 |             | Black  | King Single |
|   | S10785 | \$7,499.99 |             | Grey   | Double      |
|  | S15617 | \$7,499.99 |             | Black  | Double      |

| Image  | SKU    | Price      | Description | Colour | Sizes |
|--|--------|------------|-------------|--------|-------|
|  | S16616 | \$7,699.99 |             | Grey   | Queen |
|  | S10590 | \$7,699.99 |             | Black  | Queen |

## PRODUCT DESCRIPTION

The Icare IC333 is a premium homecare bed. While it has the features of a hospital bed, it has an appearance and practicality suited to home settings.

It provides five powered functions – backrest, leg rest, height adjustment, Trendelenburg and reverse Trendelenburg positioning – all controlled via an intuitive, backlit hand control. The bed also includes a rechargeable battery backup to maintain short-term functionality during power outages.

Its ultra-low minimum height and heavy-duty castors support safe transfers and flexible placement within the home. The breathable mattress platform, commercial-grade fabric and mattress retainer bars support comfort and hygiene over extended use.

### Who it's for

- For clients who need adjustable-height positioning and independent adjustments for safe in-home care.
- For therapists wanting a functional, durable bed to support positioning and transfers for their client.

### Key features

- **Five-function electric adjustment** – supports backrest, knee break, height adjustment, Trendelenburg and reverse Trendelenburg function
- **Ultra-low minimum height** – enhances safety, especially for those at risk of falls
- **Rechargeable battery backup** – maintains short-term operation during power outages
- **Eight heavy-duty castors** – ensures mobility and positioning flexibility
- **Backlit wipeable hand control** – easy to use day or night
- **Breathable mattress platform** – supports airflow and comfort during long-term use
- **Commercial-grade fabric** – durable and easy to maintain
- **Mattress retainer bars** – helps keep mattress securely in place

## ADDITIONAL INFORMATION

### Overall Width

Long Single - 96 cm  
King Single - 112 cm  
Double - 140 cm  
Queen - 158 cm  
Split Queen - 158 cm

### Overall Length

Long Single - 215 cm  
King Single - 215 cm  
Double - 215 cm  
Queen - 215 cm  
Split Queen - 215 cm

### Mattress Platform Width

Long Single - 90 cm  
King Single - 107 cm  
Double - 135 cm  
Queen - 152 cm  
Split Queen - 75 cm

### Mattress Platform Length

Long Single - 203 cm  
King Single - 203 cm  
Double - 203 cm  
Queen - 203 cm  
Split Queen - 203 cm

### Mattress Platform Height

Long Single - 22 cm - 66 cm  
King Single - 22 cm - 66 cm  
Double - 22 cm - 66 cm  
Queen - 22 cm - 66 cm  
Split Queen - 22 cm - 66 cm

### Overall Height

Long Single - 25 cm - 69 cm  
King Single - 25 cm - 69 cm  
Double - 25 cm - 69 cm  
Queen - 25 cm - 69 cm  
Split Queen - 25 cm - 69 cm

**Backrest Angle Adjustment**

Long Single – 0° to 62°  
King Single – 0° to 62°  
Double – 0° to 62°  
Queen – 0° to 62°  
Split Queen – 0° to 62°

**Knee Break Angle Adjustment**

Long Single – 0° to 42°  
King Single – 0° to 42°  
Double – 0° to 42°  
Queen – 0° to 62°  
Split Queen – 0° to 62°

**Trendelenburg Tilt Degree**

Long Single – 0° to 17°  
King Single – 0° to 17°  
Double – 0° to 17°  
Queen – 0° to 17°  
Split Queen – 0° to 17°

**Safe Working Load (SWL)**

Long Single -200 kg  
King Single – 200 kg  
Double -200 kg  
Queen – 0° to 17°  
Split Queen -200 kg

**Brand**

Icare Medical