



SKU:

**ICARE IC333
HOMECARE BED -
RANGE**

~~\$9,999.99~~ **\$9,999.99**

The Icare IC333 Homecare Bed combines clinical-grade functionality with a home-friendly design. Its ultra-low height, full adjustability and durable build make it a reliable option for long-term care at home.



Scan QR code to
view product online

VARIATIONS


Image	SKU	Price	Description	Colour	Sizes
	S302366	\$9,999.99		Grey	Split Queen





Image	SKU	Price	Description	Colour	Sizes
	S302365	\$9,999.99		Black	Split Queen
	S11578	\$4,299.99		Grey	Long Single
	S10589	\$4,299.99		Black	Long Single
	S11570	\$4,699.99		Grey	King Single





Image	SKU	Price	Description	Colour	Sizes
	S9907	\$4,699.99		Black	King Single
	S10785	\$7,499.99		Grey	Double
	S15617	\$7,499.99		Black	Double
	S16616	\$7,699.99		Grey	Queen

Image	SKU	Price	Description	Colour	Sizes
	S10590	\$7,699.99		Black	Queen

GALLERY IMAGES



PRODUCT DESCRIPTION

The Icare IC333 is a premium homecare bed. While it has the features of a hospital bed, it has an appearance and practicality suited to home settings.

It provides five powered functions – backrest, leg rest, height adjustment, Trendelenburg and reverse Trendelenburg positioning – all controlled via an intuitive, backlit hand control. The bed also includes a rechargeable battery backup to maintain short-term functionality during power outages.

Its ultra-low minimum height and heavy-duty castors support safe transfers and flexible placement within the home. The breathable mattress platform, commercial-grade fabric and mattress retainer bars support comfort and hygiene over extended use.

Who it's for

- For clients who need adjustable-height positioning and independent adjustments for safe in-home care.
- For therapists wanting a functional, durable bed to support positioning and transfers for their client.

Key features

- **Five-function electric adjustment** – supports backrest, knee break, height adjustment, Trendelenburg and reverse Trendelenburg function
 - **Ultra-low minimum height** – enhances safety, especially for those at risk of falls
- **Rechargeable battery backup** – maintains short-term operation during power outages
 - **Eight heavy-duty castors** – ensures mobility and positioning flexibility
 - **Backlit wipeable hand control** – easy to use day or night
- **Breathable mattress platform** – supports airflow and comfort during long-term use
 - **Commercial-grade fabric** – durable and easy to maintain
 - **Mattress retainer bars** – helps keep mattress securely in place

ADDITIONAL INFORMATION

Overall Width

Long Single - 96 cm
King Single - 112 cm
Double - 140 cm
Queen - 158 cm
Split Queen - 158 cm

Overall Length

Long Single - 215 cm
King Single - 215 cm
Double - 215 cm
Queen - 215 cm
Split Queen - 215 cm

Mattress Platform Width

Long Single - 90 cm
King Single - 107 cm
Double - 135 cm
Queen - 152 cm
Split Queen - 75 cm

Mattress Platform Length

Long Single - 203 cm
King Single - 203 cm
Double - 203 cm
Queen - 203 cm
Split Queen - 203 cm

Mattress Platform Height

Long Single - 22 cm - 66 cm
King Single - 22 cm - 66 cm
Double - 22 cm - 66 cm
Queen - 22 cm - 66 cm
Split Queen - 22 cm - 66 cm

Overall Height

Long Single - 25 cm - 69 cm
King Single - 25 cm - 69 cm
Double - 25 cm - 69 cm
Queen - 25 cm - 69 cm
Split Queen - 25 cm - 69 cm

Backrest Angle Adjustment	Long Single - 0° to 62°
	King Single - 0° to 62°
	Double - 0° to 62°
	Queen - 0° to 62°
	Split Queen - 0° to 62°
Knee Break Angle Adjustment	Long Single - 0° to 42°
	King Single - 0° to 42°
	Double - 0° to 42°
	Queen - 0° to 62°
	Split Queen - 0° to 62°
Trendelenburg Tilt Degree	Long Single - 0° to 17°
	King Single - 0° to 17°
	Double - 0° to 17°
	Queen - 0° to 17°
	Split Queen - 0° to 17°
Safe Working Load (SWL)	Long Single -200 kg
	King Single - 200 kg
	Double -200 kg
	Queen - 0° to 17°
	Split Queen -200 kg
Brand	Icare Medical

[Icare IC333 Bed User Manual](#)

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](#)

WHAT ARE THE KEY CLINICAL BENEFITS OF USING AN ADJUSTABLE BED?

HOW DO ADJUSTABLE BEDS HELP IN PREVENTING FALLS?

HOW DOES AN ADJUSTABLE BED ASSIST WITH SAFER MANUAL HANDLING FOR CARERS?

CAN ADJUSTABLE BEDS PROVIDE COMFORT AND PAIN RELIEF?

HOW DO ADJUSTABLE BEDS PROMOTE INDEPENDENCE FOR INDIVIDUALS WITH LIMITED MOBILITY?

WHAT ROLE DO ADJUSTABLE BEDS PLAY IN PRESSURE CARE?

HOW DO ADJUSTABLE BEDS FACILITATE IMPROVED BREATHING?

CAN ADJUSTABLE BEDS IMPROVE CIRCULATION?

HOW DO ADJUSTABLE BEDS HELP REDUCE SYMPTOMS OF ACID REFLUX?

WHAT IS THE TRENDELENBURG POSITION, AND WHEN IS IT USED?

WHAT IS THE REVERSE TRENDELENBURG POSITION, AND WHEN IS IT USED?

