



SKU:

ICARE IC20 MATTRESS

\$2,100.00 <u>\$2,100.00</u>

The IC20, which is in the middle range of lcare mattresses, is popular with most people. The medium feel delivers a gentle sensation while yet providing adequate support.

The IC20 medium mattress is heat and pressure responsive, which means it softens as it warms. So, when a pressure point develops, which normally becomes heated and inflamed, this portion of the mattress softens, allowing blood to flow freely again.



Scan QR code to view product online



VARIATIONS

Price **Description** Sizes Image SKU



Long Single



Long Split Queen



Long King Single



S10786 \$1,799.99

Long Double



Image SKU Price Description Sizes



Long Queen



PRODUCT DESCRIPTION

Pressure Redistribution

The IC20 mattress features the ActiveX $^{\text{m}}$ 84 material is the main comfort layer in this mattress. ActiveX $^{\text{m}}$ is a high grade of elastic type foam which is responsive to temperature around inflamed areas. ActiveX $^{\text{m}}$ will become softer as it becomes warmer. So around a pressure point it will always soften and allow the blood to circulate freely.

Solid Support Foam Base

XP Support foam is a high-density foam that is used in the bottom layer of the IC20. XP Support foam is an open cell material that breaths to allow air circulation but is high density to prevent bottoming out. This material will easily contour to any profiling bed.

Temperature Regulating Cover

The IC20 cover is made with the new Bluetec cooling fabric. This is all way stretch fabric to conform with body and mattress contour. It is under woven with a cooling yarn to regulate body temperature.

Ventilated Base Fabric

Vant fabric is a three-dimensional polypropylene material that has high air permeability and is antibacterial. We use this on the side wall of the IC20 mattress covers to greatly increase air movement.

Key features

- 200 mm Mattress Depth
- SWL Minimum 30 kg
 SWL Maximum 200 kg
- For use with the any Icare Beds
- The mattresses are also available in two more variations: soft plush IC25 and supportive firm
 IC15



ADDITIONAL INFORMATION

Sizes Long Single, Long King Single, Long Double, Long Queen, Long

Split Queen

Safe Working Load (SWL) 30 kg - 200kg

Length 203 cm

Width 90cm, 107 cm, 135 cm, 152 cm, 75 cm

Depth 20 cm

Brand Icare Medical