



SKU:

ICARE IC20 MATTRESS

~~\$2,100.00~~ **\$2,100.00**

The IC20, which is in the middle range of Icare mattresses, is popular with most people. The medium feel delivers a gentle sensation while yet providing adequate support.

The IC20 medium mattress is heat and pressure responsive, which means it softens as it warms. So, when a pressure point develops, which normally becomes heated and inflamed, this portion of the mattress softens, allowing blood to flow freely again.



Scan QR code to
view product online

VARIATIONS

SKU	Price	Description	Sizes
S10390	\$899.99		Long Single
S14399	\$2,100.00		Long Split Queen
S9909	\$999.99		Long King Single
S10786	\$1,799.99		Long Double
S10464	\$1,899.99		Long Queen

PRODUCT DESCRIPTION

Pressure Redistribution

The IC20 mattress features the ActiveX™ 84 material is the main comfort layer in this mattress. ActiveX™ is a high grade of elastic type foam which is responsive to temperature around inflamed areas. ActiveX™ will become softer as it becomes warmer. So around a pressure point it will always soften and allow the blood to circulate freely.

Solid Support Foam Base

XP Support foam is a high-density foam that is used in the bottom layer of the IC20. XP Support foam is an open cell material that breaths to allow air circulation but is high density to prevent bottoming out. This material will easily contour to any profiling bed.

Temperature Regulating Cover

The IC20 cover is made with the new Bluetec cooling fabric. This is all way stretch fabric to conform with body and mattress contour. It is under woven with a cooling yarn to regulate body temperature.

Ventilated Base Fabric

Vant fabric is a three-dimensional polypropylene material that has high air permeability and is anti-bacterial. We use this on the side wall of the IC20 mattress covers to greatly increase air movement.

Key features

- 200 mm Mattress Depth
- SWL Minimum 30 kg
SWL Maximum 200 kg
- For use with the any Icare Beds
- The mattresses are also available in two more variations: soft plush [IC25](#) and supportive firm [IC15](#)

ADDITIONAL INFORMATION

Sizes	Long Single, Long King Single, Long Double, Long Queen, Long Split Queen
Safe Working Load (SWL)	30 kg - 200kg
Length	203 cm
Width	90cm, 107 cm, 135 cm, 152 cm, 75 cm
Depth	20 cm
Brand	Icare Medical

Talk to us on the phone

Our experienced phone sales team can assist you with any product enquiries and help you make a purchase directly on the phone. We are available anytime between 9am and 5pm, Monday to Friday.

[1300 008 267](tel:1300008267)

[**WHAT ARE THE BENEFITS OF USING A PRESSURE-REDUCING MATTRESS?**](#)

[**WHAT IS THE DIFFERENCE BETWEEN A FOAM MATTRESS AND A HYBRID PRESSURE CARE MATTRESS?**](#)

[**HOW DOES AN ALTERNATING AIR MATTRESS WORK?**](#)

[**WHO CAN BENEFIT FROM USING A HYBRID PRESSURE CARE MATTRESS?**](#)

[**HOW OFTEN SHOULD FOAM PRESSURE-REDUCING MATTRESSES BE REPLACED?**](#)

[**WHAT ARE THE KEY FEATURES OF AN ALTERNATING AIR MATTRESS?**](#)

[**CAN PRESSURE CARE MATTRESSES BE USED ON ANY TYPE OF BED FRAME?**](#)

[**WHAT IS THE DIFFERENCE BETWEEN LOW AIR LOSS AND ALTERNATING AIR MATTRESSES?**](#)

[**ARE PRESSURE CARE MATTRESSES COVERED BY HEALTH INSURANCE OR NDIS?**](#)

[**HOW DO I MAINTAIN AND CLEAN AN ALTERNATING AIR MATTRESS?**](#)