

SKU:

I-CARE IC30 MATTRESS

~~\$2,199.99~~ \$2,199.99





The IC30 Gel-Infused ActiveX™ Mattress with Charcoal and Firm Edge Support is crafted to deliver exceptional comfort, support, and pressure relief for a restful night's sleep. The gel infusion regulates temperature, ensuring a cool and comfortable sleeping surface, while the charcoal absorbs odours and moisture, promoting a clean and hygienic environment.

ActiveX™ technology is temperature- and pressure-sensitive, offering an ideal balance of immersion and envelopment. This results in a sleep surface that provides outstanding support, comfort, and high-pressure care.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S325734	\$1,199.99		Long Single
	S314120	\$1,299.99		King Single
	S322409	\$2,099.99		Double
	S325755	\$2,199.99		Queen

PRODUCT DESCRIPTION

- Pressure-relieving material
- Flip Free Design
- Vant ventilated sides and base fabric with Bluetec temperature-regulating cover
- Machine washable covers
- User Weight: 30 kg - 200 kg
- Waterlow ≤ 20

ADDITIONAL INFORMATION

Brand	I-Care Medical
Size	Double, King Single, Queen, Long Single
Measurements	Double: 1900 x 1530
	King Single: 2030 x 1070
	Queen: 2030 x 1520
	Long Single: 2030 x 900