



SKU:

DRIVE MEDICAL INDIANA LIFT CHAIR (SINGLE MOTOR)

\$2,299.99

The Drive Medical Indiana Lift Chair (Single Motor) provides powered assistance for sitting and standing. With plush fabric and simple controls, it delivers supportive comfort for everyday use. Available in Petite and Standard sizes.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colours	Sizes
	S11881	\$2,299.99		Mushroom	Standard
	S15398	\$2,299.99		Mushroom	Petite
	S12781	\$2,299.99		Grey Graphite	Standard
	S15397	\$2,299.99		Grey Graphite	Petite

PRODUCT DESCRIPTION

The Drive Medical Indiana Lift Chair (Single Motor) is a compact and practical option for individuals who need help transitioning between sitting and standing positions. Its single motor mechanism allows for smooth, quiet movement, operated using a simple two-button handset. The chair gently reclines for relaxation and rises to support safe standing.

Upholstered in durable 'Sofia' fabric, available in Graphite Grey or Mushroom, the Indiana combines modern styling with everyday functionality.

The ergonomically contoured seat, fully padded chaise and supportive backrest enhance comfort during extended sitting.

A side pocket provides a convenient place to keep essentials within reach.

Who it's for

- For clients who need gentle assistance with standing and reclining in daily use.
- For therapists wanting a reliable lift chair with simple operation and ergonomic design for their client.

Key features

- **Single motor lift and recline** – supports safe transitions from sitting to standing
- **Two-button handset** – easy to use with minimal effort
- **Contoured seat and backrest** – promotes posture and comfort
- **Fully padded chaise** – enhances leg support and relaxation
- **Sofia fabric upholstery** – soft, durable and available in two modern colours
- **Discreet and quiet operation** – ideal for shared or home environments
- **Side storage pocket** – keeps personal items close at hand
- **Fire retardant materials** – complies with safety standards
- **Available in Petite and Standard sizes** – offers a better fit for different needs

ADDITIONAL INFORMATION

Brand	Drive Medical
Max SWL	150KG
Width	81cm (32") ST 77cm (30") PE
Depth	96cm (37") ST 88cm (34") PE
Seat Width	51cm (20") ST 47cm (18") PE
Seat Height	48cm (19") ST 44cm (17") PE
Seat Depth	53cm (21") ST 48cm (19") PE
Colours	Mushroom, Grey Graphite
Sizes	Petite, Standard