



SKU:

DRIVE MEDICAL INDIANA RECLINER LIFT CHAIR


~~\$1,999.99~~ \$1,999.99

The Indiana posture care recliner lift chair is both sophisticated and comfortable, providing outstanding value for your money. The lift chair is available in standard or petite sizing options for your comfort.

The Indiana Rise & Recliner is furnished in plush Sophia fabric and is available in Graphite, Russet and Mushroom color options.

Product features include

- New luxurious Sofia fabric.
- 2 button handset with safety feature.
- Reclines to snooze position.
- Convenient side pocket.
- The footrest can be elevated.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colours	Sizes
	S11881	\$1,899.99		Mushroom	Standard
	S15398	\$1,999.99		Mushroom	Petite
	S12781	\$1,999.99		Grey Graphite	Standard
	S15397	\$1,999.99		Grey Graphite	Petite

PRODUCT DESCRIPTION

The Indiana posture care recliner lift chair is stylish and comfy, and it offers excellent value for money. For your convenience, the reclining sofa is available in regular or petite sizes, and it would make a terrific lounge in any living area.

The Indiana Recliner Lift Chair is upholstered in soft Sophia fabric and comes in three color options: Graphite, Russet, and Mushroom.

The Indiana Rise lounge chair has a full length chaise. It has a soft and comfortable seat and back cushion for added comfort. The lift chair will rise and tilt you to your feet. The warranty periods is 2 years and the weight capacity is 150kg.

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

Brand	Drive Medical
Max SWL	150KG
Width	81cm (32") ST 77cm (30") PE
Depth	96cm (37") ST 88cm (34") PE
Seat Width	51cm (20") ST 47cm (18") PE
Seat Height	48cm (19") ST 44cm (17") PE
Seat Depth	53cm (21") ST 48cm (19") PE
Colours	Mushroom, Grey Graphite
Sizes	Standard, Petite