



**SKU:**

## **DRIVE MEDICAL INDIANA LIFT CHAIR (SINGLE MOTOR)**

**\$2,299.99**

The Drive Medical Indiana Lift Chair (Single Motor) provides powered assistance for sitting and standing. With plush fabric and simple controls, it delivers supportive comfort for everyday use. Available in Petite and Standard sizes.



Scan QR code to  
view product online

## VARIATIONS

Image	SKU	Price	Description	Sizes	Colours
	S11881	\$2,299.99		Standard	Mushroom
	S15398	\$2,299.99		Petite	Mushroom
	S12781	\$2,299.99		Standard	Grey Graphite
	S15397	\$2,299.99		Petite	Grey Graphite

## PRODUCT DESCRIPTION

The Drive Medical Indiana Lift Chair (Single Motor) is a compact and practical option for individuals who need help transitioning between sitting and standing positions. Its single motor mechanism allows for smooth, quiet movement, operated using a simple two-button handset. The chair gently reclines for relaxation and rises to support safe standing.

Upholstered in durable 'Sofia' fabric, available in Graphite Grey or Mushroom, the Indiana combines modern styling with everyday functionality.

The ergonomically contoured seat, fully padded chaise and supportive backrest enhance comfort during extended sitting.

A side pocket provides a convenient place to keep essentials within reach.

### Who it's for

- For clients who need gentle assistance with standing and reclining in daily use.
- For therapists wanting a reliable lift chair with simple operation and ergonomic design for their client.

### Key features

- **Single motor lift and recline** – supports safe transitions from sitting to standing
- **Two-button handset** – easy to use with minimal effort
- **Contoured seat and backrest** – promotes posture and comfort
- **Fully padded chaise** – enhances leg support and relaxation
- **Sofia fabric upholstery** – soft, durable and available in two modern colours
- **Discreet and quiet operation** – ideal for shared or home environments
- **Side storage pocket** – keeps personal items close at hand
- **Fire retardant materials** – complies with safety standards
- **Available in Petite and Standard sizes** – offers a better fit for different needs

## ADDITIONAL INFORMATION

<b>Safe Working Load (SWL)</b>	150 kg ST
	150 kg PE
<b>Overall Width</b>	77 cm ST
	81 cm PE
<b>Overall Depth</b>	96cm ST
	88cm PE
<b>Overall Height</b>	108 cm ST
	103 cm PE
<b>Seat Width</b>	51 cm ST
	47 cm PE
<b>Seat Depth</b>	53 cm ST
	48 cm PE
<b>Seat to Floor Height</b>	48 cm ST
	44 cm PE
<b>Number of Motors</b>	Twin
<b>Sizes</b>	Petite, Standard
<b>Colours</b>	Mushroom, Grey Graphite
<b>Brand</b>	Drive Medical