



SKU:

# INOGEN ROVE 6 PORTABLE OXYGEN CONCENTRATOR

**\$4,699.99 \$4,699.99** 

The Inogen Rove 6 Portable Oxygen Concentrator provides dependable oxygen therapy in a compact design. It supports daily activities and travel with a lightweight build, six flow settings and quiet operation, enhancing comfort and mobility.



Scan QR code to view product online



# **VARIATIONS**

Image	SKU	Price	Description	Type of battery
	S330512 \$	4,499.99	8 cell	
	S330509 \$	4,699.99	16 cell	



### PRODUCT DESCRIPTION

The Inogen Rove 6 Portable Oxygen Concentrator is designed to support those requiring supplemental oxygen therapy during daily living and travel.

Weighing only 2.2 kg with a single battery, the Inogen Rove 6 is easy to carry and transport. It uses pulse-dose delivery to synchronise oxygen output with the person's breathing pattern, ensuring efficient therapy during activity or rest.

The concentrator offers six adjustable oxygen flow settings to cater to a range of respiratory needs. The bright digital display ensures intuitive operation, while the device's compact design allows for use at home, in clinical settings, or while travelling.

It operates on rechargeable battery or AC/DC power, offering flexibility and convenience across different environments.

Medical-grade oxygen is delivered through a nasal cannula, with Intelligent Delivery Technology ensuring precise, rapid breath detection and responsive therapy.

The Inogen Rove 6 is supplied with AC power supply, DC power cable, carry bag and user manual.

### Who it's for

- For clients requiring supplemental oxygen therapy seeking a portable, lightweight solution for daily activities and travel.
- For therapists prescribing oxygen therapy who need a reliable, user-friendly device for their clients.

## **Key features**

- **Compact and lightweight design** weighs only 2.2 kg with a single battery, supporting easy portability
- Versatile use suitable for use at home, in clinical settings, and during transport
- Adjustable oxygen flow settings provides six high-flow options to suit individual respiratory needs
- **Pulse dose delivery** synchronises oxygen delivery with breathing pattern for efficient therapy
- Bright digital display easy-to-read interface for intuitive operation
- **Multiple power options** operates on rechargeable battery or standard AC/DC power, offering flexibility across environments
- Medical-grade oxygen delivery uses a nasal cannula to supply oxygen-rich air by removing nitrogen



### **ADDITIONAL INFORMATION**

Weight 2.2 kg (single battery)

2.6 kg (with extended-life battery)

**Length** 18.3 cm

Width 8.3 cm

Height 20.5 cm

Oxygen Flow Pulse dose

**Power** AC Power Supply 100-240V

**Battery** Single battery up to 6.25 hours

Double battery up to 12.75 hours

Noise level 39 dBA\*

**Operation** Simple control functions and easy-to-read LCD display.

Intended use Designed for 24/7 use and conforms to all applicable FAA

requirements for carriage and use on board aircraft

**Brand** Inogen

Oxygen Flow Pulse Dose