




SKU:

## **INOGEN ROVE 6 PORTABLE OXYGEN CONCENTRATOR**

**~~\$6,199.99~~ \$6,199.99**

The Inogen Rove 6 is an ideal solution for patients looking for a Portable Oxygen Concentrator that maximizes their range and flexibility for oxygen delivery in a system that is quiet, easy to use and easy to maintain. The Inogen Rove 6 is designed to empower your patients with the peace of mind that comes with feeling ready to continue their normal activities - like spending quality time with family and friends, enjoying the great outdoors and traveling. With 6 flow settings, the Inogen Rove 6 is the most flexible POC we've ever offered, delivering up to 1,260 mL of medical-grade oxygen per minute.



Scan QR code to  
view product online

## VARIATIONS

**Image**

**SKU**

**Price**

**Description**

**Type of battery**



S330512 \$5,899.99

8 cell



S330509 \$6,199.99

16 cell

## PRODUCT DESCRIPTION

The Inogen Rove 6 concentrator weighs just 2.2kg (single battery). It can be charged at home or on the go—giving you the freedom and independence to where you want, when you want, and at a moment's notice! The Inogen Rove 6 system was designed to be lightweight, quiet and minimise service requirements.

### Key Features

- Compact and lightweight unit – it weighs just 2.2 kg with single battery
- Used in the home, institution, and transport modalities
- High flow oxygen settings from 1-6
- Bright easy-to-read digital display
- Delivers oxygen via Pulse Dose mechanisms. Pulse Dose technology utilises an oxygen conserver to deliver oxygen based on inhalation rates, allowing the proper amount of oxygen to be delivered during all phases of daily activity and during rest
- Powered by a rechargeable battery or standard AC or DC power source
- Provide a high concentration of supplemental oxygen to patients requiring respiratory therapy on a prescriptive basis
- Deliver medical grade oxygen via a nasal cannula, removing nitrogen to deliver oxygen-rich air
- Intelligent Delivery Technology found in all Inogen products utilises triggering sensitivity to quickly detect a breath and deliver oxygen

### The Inogen Rove 6 comes with the following:

- AC power supply
- DC power cable
- Carry Bag
- User Manual

## ADDITIONAL INFORMATION

<b>Brand</b>	Inogen
<b>Weight</b>	2,2 kg. (4.8 lbs.) (includes single battery) 2,6 kg. (5.8 lbs.) (with extended-life battery)
<b>Length</b>	18.3 cm
<b>Width</b>	8.3 cm
<b>Height</b>	20.5 cm
<b>Oxygen flow</b>	Pulse dose delivery system Six flow settings: 1-6
<b>Power</b>	AC Power Supply 100-240V, 50-60Hz (auto sensing 2.0 – 1.0A to allow worldwide use) DC Power Cable: for mobile use in a vehicle
<b>Battery</b>	Single Battery up to 6.25 hours Double Battery up to 12.75 hours
<b>Noise level</b>	39 dBA*
<b>Operation</b>	Simple control functions and easy-to-read LCD display.
<b>Intended use</b>	Designed for 24/7 use and conforms to all applicable FAA requirements for carriage and use on board aircraft