



INOGEN ONE G5 PORTABLE OXYGEN CONCENTRATOR

\$4,649.99

Inogen innovation sets a new standard in oxygen delivery with the G5 portable oxygen concentrator. The Inogen One G5 offers the most oxygen per kilogram for a portable oxygen concentrator. With an efficient and lightweight design, the Inogen One G5 offers fluctuating flow settings from 1-6 and is engineered for 24-hour use, boosting your independence.

It can be charged at home or on the go giving you the freedom of oxygen. Anytime. Anywhere. The Inogen One G5 is your one solution, your single solution, for oxygen at home or away.



SKU: S19446

Categories: [Oxygen Concentrator](#), [Portable Oxygen Concentrators](#)

GALLERY IMAGES



PRODUCT DESCRIPTION

The Inogen One G5 portable oxygen concentrator weighs just 2.2kg with its battery, making it both smaller and lighter than the Inogen One G3. It can be charged at home or on the go—giving you the freedom and independence to where you want, when you want, and at a moment's notice!

Key features include:

- Compact and lightweight unit – it weighs just 2.2kg with single battery
 - High flow oxygen settings from 1-6
 - Intelligent pulse oxygen delivery
- Concentrated oxygen for 24 hour use
- Easy-to-read symbolized LCD display
 - Simple control functions
- User-replaceable batteries and sieve beds
 - Warranty period of 3 years

The Inogen One G5 comes with the following:

1. 8 cell single battery
2. AC power supply
3. DC power cable
4. Custom carrying case
5. Shoulder strap

6. Tubing
7. Filter

To compare the Inogen G5 with the G3, read our [blog](#).

Please note that a prescription is required for you to purchase medical oxygen.

ADDITIONAL INFORMATION

| | |
|--------------------|---|
| Brand | Inogen |
| Weight | 2.2kg (4.7 lbs) includes single battery |
| Length | 18.26cm (7.19 inches) |
| Width | 8.28cm (3.26 inches) |
| Height | 20.7cm (8.15 inches) |
| Oxygen flow | Pulse dose delivery system Six flow settings: 1-6 |
| Power | AC power supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC power cable: for mobile use in car |
| Battery | Duration (single battery): Up to 6.5 hours. Recharge up to 3 hours with AC or DC power. Duration (double battery): Up to 13 hours. Recharge up to 6 hours with AC or DC power. |
| Noise level | 38 dBA* |
| Operation | Simple control functions and easy-to-read LCD display Inogen Connect App |
| Use | Designed for 24/7 use and conforms to all applicable FAA requirements for POC use on board an aircraft. |