

SKU: S15324




INOGEN ONE G4 OXYGEN CONCENTRATOR

~~\$6,199.99~~ **\$6,199.99**

Introducing the next innovation to portable oxygen delivery: the Inogen One G4, capable of full 24/7 oxygen delivery in an unbelievably small package!.

The Inogen One G4 is officially the smallest and lightest portable oxygen concentrator ever released by Inogen. Weighing in at 1.27kg with the single battery, the Inogen One G4 is more than 40% lighter than the Inogen One G3 and it's about half the size! It even weighs less than some of the most commonly used portable oxygen tanks.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Type of battery
	S15324	\$5,899.99		4 cell
	S305270	\$6,199.99		8 cell

PRODUCT DESCRIPTION

For such a small unit, the Inogen One G4 is incredibly powerful. It's a pulse flow delivery system, which means it does not offer continuous flow, but it offers settings 1-3. It is equipped with Inogen's standard Intelligent Delivery Technology to ensure the unit detects your breath to delivery your oxygen with incredible accuracy. The G4 is the smallest and lightest unit in its class, and you don't sacrifice functionality or battery life for the small size.

- Unbelievably Compact & Lightweight
- Ultra-Quiet Operation
- Intelligent Pulse Oxygen Delivery
- Designed for 24 Hour Use

What Comes With The Inogen One G4?

- Inogen One G4 Portable Oxygen Concentrator
- Inogen 4 cell battery
- AC Power Supply
- DC Power Supply
- Custom Carrying Case
- Shoulder Strap
- Tubing
- Filter
- Manual
- 3-Year Warranty

Doctor's prescription required for dispatch of any oxygen concentrator.

ADDITIONAL INFORMATION

Weight	1.27kg
Size	Length: 15cm X Width: 6.8cm X Height: 18cm (includes single battery)
Flow Rate	Three flow settings 1, 2, 3 – Pulse dose delivery system
Power	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
Battery	Duration (single battery): Up to 3 hours Inogen G4 4 Cell Battery
Noise Level	40 dBA (at flow setting 2)
Warranty	3 year warranty
Operation	Simple control functions
Use	Designed for 24/7 use