




SKU:

INVACARE COLIBRI MOBILITY SCOOTER

\$1,999.99


The Invacare Colibri scooter is a compact and lightweight mobility scooter that is ideal for individuals who enjoy an independent lifestyle. It is a small-sized scooter that can be easily maneuvered and stored in small spaces. This scooter has a unique 'LiteLock' system to allow a simple, tool-free disassembly. The Invacare Colibri is designed with a user-friendly interface, making it easy to operate. It has a sleek and modern design, with a range of customisable features to suit individual needs and preferences. Users can select from a range of trendy colour panels to match individual style, swapping them according to your needs or when a refresh is needed.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colours
	B10007	\$1,999.99		Orange
	B10009	\$1,999.99		White
	B10008	\$1,999.99		Jasper Red
	B10005	\$1,999.99		Sapphire Blue

Image	SKU	Price	Description	Colours
	B10006	\$1,999.99		Jade Green

PRODUCT DESCRIPTION

The Invacare Colibri scooter features a sturdy and durable design, with a weight capacity of up to 136 kg. It comes with a comfortable seat, which has multiple adjustments to suit user needs and provide maximum comfort. The Colibri has a range of up to 16 km on a single charge, making it ideal for short trips and errands. It has a top speed of 8 km/h, and its front and rear suspension system ensures a smooth and comfortable ride.

Key Features

- Compact and lightweight design
- Invacare LiteLock system ensures the scooter can easily be taken apart without the requirement of tools
- User-friendly interface
- Comfortable seat and armrests
- Range of up to 16 km- Please note this is influenced by external factors
- Top speed of 8 km/h
- Customizable features to suit individual needs and preferences
- Invacare Scooter range comes with FREE 24/7 road assistance
- Multiple colours to choose from which match your style

ADDITIONAL INFORMATION

Brand	Invacare
Model	Colibri
Weight Capacity / Max SWL	136 kg
Total Length	101 cm
Total Width	50.5 cm
Total Height	84 cm
Floor to Seat Height	37 cm – 42 cm
Backrest Height	38.5 cm
Backrest Width	46.5 cm
Seat Depth	40 cm
Front/Rear Tyre Type	Solid
Front/Rear Tyre Size	7.8' inch
Ground Clearance	5 cm
Turning Diameter	220 cm
Max Safe Slope	6 degrees
Battery Capacity/Size	2 x 12V/18Ah
Motor Capacity	200 W
Charger Type	2AH
Maximum Driving Speed	8Km/h
Maximum Driving Distance/Range	Up to 16Km

Carton Size (W x L x H)

61 cm x 116 cm x 71 cm