



**SKU:**

## **INVACARE COLIBRI MOBILITY SCOOTER**

**\$1,999.99**

The Invacare Colibri Mobility Scooter is a lightweight, compact scooter designed for independence and everyday mobility. Its unique LiteLock system enables easy tool-free disassembly, so it's easy to store and transport.



Scan QR code to  
view product online

**VARIATIONS**






Image	SKU	Price	Description	Colours
	B10007	\$1,999.99		Orange
	B10009	\$1,999.99		White
	B10008	\$1,999.99		Red
	B10005	\$1,999.99		Blue

Image	SKU	Price	Description	Colours
	B10006	\$1,999.99		Green

## PRODUCT DESCRIPTION

The Invacare Colibri Mobility Scooter offers a practical, user-friendly solution for those seeking convenient short-distance travel. With a modern design and intuitive controls, it suits both indoor and outdoor environments.

This small-sized scooter is easy to steer and manoeuvre in tight spaces. It features a soft padded swivel seat with width-adjustable armrests, and offers multiple seating adjustments for added comfort.

With a range of up to 16 km on a full charge and a top speed of 8 km/h, the Colibri supports short trips with ease.

The innovative LiteLock system enables quick, tool-free disassembly into manageable parts, ideal for car transport or storage. Colour panels can be swapped out to suit individual preference, offering a fresh look whenever desired.

All Invacare scooters include 12 months of complimentary 24/7 roadside assistance upon registration with the manufacturer.

### Who it's for

- For clients wanting a compact scooter that's easy to dismantle and store.
- For therapists recommending a lightweight scooter with adjustable seating and personalisation options.

### Key features

- **LiteLock tool-free disassembly** – supports easy transport or storage
- **Compact, lightweight frame** – for easier manoeuvrability in small spaces
- **Padded swivel seat with adjustable armrests** – promotes comfort and safe transfers
- **Range up to 16 km with 8 km/h speed** – ideal for short trips and errands
- **Colour panel customisation** – allows personal style preferences
- **User-friendly controls** – designed for intuitive operation
- **12-month 24/7 roadside assistance** – free with registration through Invacare

## ADDITIONAL INFORMATION

**Brand**

Invacare