




SKU:

INVACARE COMET ULTRA MOBILITY SCOOTER

\$6,299.99

The Invacare Comet Ultra Mobility Scooter offers a secure, supportive, and powerful driving experience without compromising on style or functionality. The Comet Ultra will get you to your destination quickly and efficiently, in complete comfort and control. The Comet Ultra is designed with improved safety and stability with a safe working load up to 220Kg. Developed with the highest quality components and safety features, the Comet Ultra makes a trustworthy companion no matter where you go. Furthermore, its flexibility and adaptability means you can customise your scooter to suit your unique style and individual preference with a variety of shroud colours.



Scan QR code to
view product online

VARIATIONS








Image	SKU	Price	Description	Colours
	B10019	\$6,299.99		Silver
	B10017	\$6,299.99		Green
	B10018	\$6,299.99		Sapphire Blue
	B10048	\$6,299.99		Electric Blue

Image	SKU	Price	Description	Colours
	B10020	\$6,299.99		Red
	B10021	\$6,299.99		Crystal
	B10016	\$6,299.99		White

PRODUCT DESCRIPTION

The Invacare Ultra Heavy-Duty scooter boasts a robust design that can handle a weight capacity of up to 220kg. The scooter is equipped with full suspension, which provides a smooth ride and increased comfort, even on bumpy surfaces. The special contoured 'ultra-seat' has been developed in conjunction with therapists to ensure an optimum, ergonomic seating position and increased comfort. Multiple adjustments to the seat and steering column to suit user needs and desired comfort.

Key Features

- **Weight capacity:** The Invacare Comet Ultra scooter has a weight capacity of up to 220kg
- **Safety:** The scooter comes with various safety features, including a headlight, taillight, and turn signals, ensuring that users remain visible and safe when traveling
- **Powerful motor:** The Comet Ultra scooter has a powerful motor that can reach speeds of up to 10Km/h and a maximum driving distance of 49Km
- **Comfortable seat:** The scooter features a comfortable and ergonomic seat that can be adjusted for height, depth, and angle, ensuring optimal comfort and support for the user.
- **Adjustable Steering Column:** The scooter's steering column/controls can be adjusted to accommodate different arm lengths and positions, ensuring comfortable and easy steering.
- **Large wheels:** The Comet Ultra scooter features large, air-filled tires that provide enhanced traction and stability on different terrains, making it suitable for both indoor and outdoor use.
- **Durability:** The Comet Ultra Heavy Duty is built with high-quality materials, including a robust steel frame, which provides superior durability and longevity.
- **Colours:** Multiple shroud colour options

ADDITIONAL INFORMATION

Brand	Invacare
Model	Comet Ultra HD
Weight Capacity / Max SWL	220 kg (485 lbs)
Empty Weight	148 kg
Total Length	145 cm
Driving Unit Width	66.5 cm - 68.5 cm
Total Height	129 cm
Floor to Seat Height	44 cm - 51 cm
Backrest Height	51.5 cm
Backrest Width	66 cm
Seat Depth	53.5 cm
Front/Rear Tyre Type	Pneumatic
Front/Rear Tyre Size	13' inch
Ground Clearance	10 cm
Turning Diameter	275 cm
Max Safe Slope	10 degrees
Battery Capacity/Size	2 x 12V/75Ah
Motor Capacity	650W/1800W
Charger Type	8AH
Maximum Driving Speed	10Km/h

**Maximum Driving
Distance/Range**

49Km

Carton Size (W x L x H)

73 cm x 174 cm x 74 cm