



INVACARE ISA STAND ASSIST AID XPLUS

~~\$4,595.99~~ \$3,995.99


Our new versatile Invacare Stand Assist lifter has been designed to suit a range of users who require assistance during a transfer but can weight bear. The new Invacare Stand Assist lifter provides safety, comfort and a more pleasant experience for clients. ISA is the product for those independent individuals who are in need of some assistance.

Suitable for people who have some weight bearing ability, but are unable to raise themselves unaided. Ergonomic Hand Grips enable positioning of the hands in a natural, comfortable and intuitive way.

Extendable Lifting Arm adapts the lifter to the person's height, condition and preference.

Lower Leg Support maintains the lower legs with maximum support.


Can be assembled and dismantled without tools to achieve ease of transportation.



SKU: S20189

Categories: [Floor Hoist and Lifters](#), [Patient Transfer Equipment](#)

GALLERY IMAGES



PRODUCT DESCRIPTION

As a caregiver, usability is a fundamental element, which is why the new ISA has a reduced, curved

overall base dimension and reduced product weight. In many cases, the caregiver may be a family member, and therefore the ease of use is imperative. Dismantling this product can be done without any tools, and thus transportation can be achieved fluidly.

ADDITIONAL INFORMATION

Base Length	1280mm
Base Width	640mm
Inside Base width	1020mm
Turning Diameter	1260mm
Machine Weight	52kg
Max SWL	200 kg
Brand	Invacare
Available on Queensland Medical Aids Subsidy Scheme	Yes
Ground Clearance	35mm
Lifts per Charge	82
Lifting Range	505/675mm