




SKU:

INVACARE LEO MOBILITY SCOOTER

\$3,599.99

The Invacare Leo is a safe, comfortable ride with a solid trusted design and quality manoeuvrability. The Invacare Leo is an optimal 4-wheeled scooter with a stylish and sporty finish. Being super compact, the Leo is a great option for scooting around the shops and can handle the tight corners around the home. Safety is a key feature with functions such as an intelligent lighting system. Available in two colours, discover the high standards of design and safety engineering with the Leo mobility scooter. The Invacare mobility scooter range comes with FREE 24/7 road assistance for the first 12 months through Invacare. Apply by registering [here](#).



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Colour
	S11107	\$3,599.99		Red
	S12768	\$3,599.99		Blue

PRODUCT DESCRIPTION

The Leo is a compact, safe, and stylish scooter that provides ultimate comfort and manoeuvrability.

Key features include:

- Intuitive control panel that is easy to operate – The control panel has an accurate easy-to-read battery display level, speed dial, horn button, hazard light button, left and right indicators, headlight button, and thumb levers for forwards/backwards
- Intelligent lighting system includes “Xenon-style” lights and indicators that automatically switch off. The scooter has rear brake lights to warn others the brake is engaged/slowing the scooter down
- 2-step disengaging lever – This unique lever prevents the scooter accidentally freewheeling
- Strong motor ensures reliable powerful performance, with a 38Km maximum driving distance, max speed of 8km/h and a safe working load of 136kg
- Luxury swivel seat – The soft padding in the seat and the backrest makes sitting for lengths of time a lot more comfortable. The white stitch work through the seat also provides a stylish edge to the overall look of the scooter. Fully adjustable for extra comfort, the seat also provides improved access to get on and off the scooter
- Flip up arm rests makes getting in and out of the seat easier
- Foot pads on front guards provides ample space on floor area for extra leg room and storage
- Simple charging system

- Easy to transport and dismantle – Longer trips or holidays are hassle free as you can quickly assemble and dismantle the Leo without using tools. It can fit in most car boots
- Easy-to-access battery and the electronics – You simply remove just one part/shell casing making it more convenient for maintenance checks and servicing
- Front carry basket for small items

ADDITIONAL INFORMATION

Brand	Invacare
Model	Leo
Weight Capacity / Max SWL	136kg (299 lbs)
Empty Weight	59 (exc. Batteries) & 83Kg (inc. Batteries)
Total Length	122 cm
Driving Unit Width	59 cm
Total Height	98 cm - 123 cm
Floor to Seat Height	46 cm -56 cm
Backrest Height	47.5 cm
Backrest Width	47 cm
Seat Depth	41 cm
Front/Rear Tyre Type	Pneumatic
Front/Rear Tyre Size	10' inch
Ground Clearance	7 cm
Turning Diameter	262 cm
Max Safe Slope	10 degrees
Battery Capacity/Size	2 x 12V/36Ah
Motor Capacity	240W
Charger Type	5AH
Maximum Driving Speed	8 km/h

**Maximum Driving
Distance/Range**

38 km

Carton Size (W x L x H)

66 cm x 143 cm x 86 cm