


SKU:

JAY BASIC CUSHION

~~\$249.99~~ **\$249.99**



The Jay Basic Cushion provides durable comfort and mild postural support for people at low risk of skin breakdown. With its contoured moulded foam base and moisture-resistant cover, it's an ideal choice for daily use.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Sizes
	S17636	\$209.99		16" x 16"
	S2946	\$209.99		16" x 18"
	S1767	\$249.99		18" x 16"
	S2945	\$249.99		18" x 18"

Image	SKU	Price	Description	Sizes
	S2947	\$249.99		20" x 18"
	S7863	\$249.99		20" x 20"

PRODUCT DESCRIPTION

The Jay Basic Cushion features a contoured, moulded foam base that provides mild lateral and posterior stability for individuals at low risk of pressure injury. Its bevelled base design accommodates seat slings, promoting a stable seating surface for those with intact skin integrity and independent weight-shifting ability.

Constructed from durable foam, the cushion is lightweight and easy to handle.

The moisture-resistant cover protects the foam from fluid build-up, while the non-slip bottom helps keep the cushion securely in place during use and transfers.

Who it's for

- For clients at low risk of pressure injury who need a lightweight, comfortable cushion for daily use.
- For therapists prescribing a practical cushion that offers mild postural support and prioritises skin health.

Key features

- **Moulded foam base with mild contouring** – provides comfort and gentle lateral and posterior stability
- **Bevelled base design** – accommodates seat slings for a stable seating surface
- **Moisture-resistant cover with non-slip bottom** – protects foam from fluid build-up and keeps cushion in position
- **Lightweight construction** – facilitates easy handling and transport

ADDITIONAL INFORMATION

Height	2"
Brand	Jay
Manufacturer's part number	J300