




SKU:

JAY FUSION CUSHION

~~\$1,339.99~~ \$1,339.99

The Jay Fusion Cushion offers effective skin protection and positioning for people at high risk of pressure injury. Its multi-layered contoured foam base and Jay Flow fluid inserts conform to each person's shape, providing reliable pressure redistribution, comfort and stability.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Sizes
	S9810	\$1,099.99		16" x 16"
	S12144	\$1,099.99		16" x 18"
	S15718	\$1,199.99		16" x 20"
	S9807	\$1,199.99		18" x 16"




Image	SKU	Price	Description	Sizes
	S9808	\$1,199.99		18" x 18"
	S11488	\$1,199.99		18" x 20"
	S7718	\$1,099.99		20" x 16"
	S11231	\$1,099.99		20" x 20"





Image	SKU	Price	Description	Sizes
	S307087	\$1,099.99		20" x 22"
	S14794	\$1,199.99		22" x 18"
	S13499	\$1,339.99		22" x 20"
	S331446	\$1,199.99		22" x 22"




Image	SKU	Price	Description	Sizes
	S303594	\$1,199.99		24" x 18"
	S303246	\$1,199.99		24" x 20"
	S9809	\$1,099.99		20" x 18"
	S15718	\$999.99		16x20"





Image	SKU	Price	Description	Sizes
	S11900	\$1,099.99		17x17"
	S12011	\$1,099.99		17x18"
	S10985	\$1,099.99		17x19"
	S12413	\$1,099.99		17x21"



Image	SKU	Price	Description	Sizes
	S9807	\$1,149.99		18x16"
	S12280	\$1,149.99		18x17"
	S9808	\$1,149.99		18x18"
	S11508	\$1,149.99		18x19"




Image	SKU	Price	Description	Sizes
	S11488	\$1,149.99		18x20"
	S12381	\$1,149.99		19x18"
	S10955	\$1,149.99		19x19"
	S11287	\$1,149.99		19x21"




Image	SKU	Price	Description	Sizes
	S7718	\$1,149.99		20x16"
	S9809	\$1,149.99		20x18"
	S11231	\$1,149.99		20x20"
	S11767	\$1,149.99		20x21"

Image	SKU	Price	Description	Sizes
	S307087	\$1,149.99		20x22"

PRODUCT DESCRIPTION

The Jay Fusion Cushion is designed for individuals who need advanced skin protection and postural support during prolonged sitting. Its multi-layered, contoured closed-cell foam base features a soft top layer that enhances positioning and comfort, while the femoral transition and rear wall support the pelvis and trochanters to help maintain fluid placement and stability.

The pelvic loading area is proportioned using average skeletal measurements, promoting even immersion and forward and lateral stability across a range of body shapes. The Jay Flow fluid inserts adjust to individual shapes and provide consistent positioning beneath bony prominences, offering long-term skin protection with minimal maintenance.

The dual-cover system includes a stretch, antimicrobial outer cover and a fluid-resistant inner cover with Aqua-guard zip and anti-wicking thread to protect the foam, supporting hygiene and durability in everyday use.

Who it's for

- For clients at high risk of pressure injury who need reliable comfort and postural support during extended sitting.
- For therapists seeking an adjustable cushion that balances pressure management, stability and ease of maintenance for clients with complex seating needs.

Key features

- **Multi-layered, contoured closed-cell foam base** – offers adjustable support and helps manage fluid migration
- **Soft top layer with femoral transition and rear wall** – enhances postural alignment and supports the pelvis and trochanters
- **Pelvic loading area based on skeletal averages** – promotes immersion and stability across different body shapes
- **Jay Flow fluid inserts** – conform to each person's shape and maintain consistent positioning beneath bony prominences
- **Dual-cover system** – stretch outer cover with antimicrobial properties and water-resistant inner cover for hygiene and durability

ADDITIONAL INFORMATION

Width	14" - 24"
Depth	14" - 24"
Product Weight Max	2.4 kg (16" x 18")
Cover Type	X-static® silver thread stretch. Incontinent & microclimatic options. Standard with 2 covers.
SWL	100kg to 149kg