

SKU:

JAY J2 CUSHION

~~\$1,349.99~~ \$1,349.99



The Jay J2 Cushion provides advanced skin protection and postural support for people at high risk of pressure injury. With a pre-contoured foam base and Jay Flow fluid tripad, it delivers reliable pressure redistribution and stability during prolonged sitting.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Sizes
	S2940	\$1,099.99		16" x 16"
	S8200	\$1,099.99		16" x 18"
	S2280	\$1,099.99		18" x 16"
	S1513	\$1,099.99		18" x 18"

Image	SKU	Price	Description	Sizes
	S2926	\$1,099.99		18" x 20"
	S318714	\$1,349.99		20" x 18"

PRODUCT DESCRIPTION

The Jay J2 Cushion is a pre-contoured foam cushion designed for individuals at high risk of pressure injury who need enhanced comfort and postural support. Its firm foam base includes an anatomical well that promotes stable positioning and effective pressure redistribution.

The cushion's Jay Flow fluid tripad conforms to the user's shape, adjusting to different sitting positions to help protect the skin from breakdown by redistributing pressure beneath bony prominences. The fluid tripad provides up to 3" of loading for effective immersion and envelopment. Constructed from closed-cell foam, the base is easily modified to accommodate each person's unique needs.

The Jay J2 Cushion also includes an Air Exchange cover, which is breathable and moisture-wicking, with optional accessories available for incontinence protection and solid seat compatibility.

Who it's for

- For clients at high risk of pressure injury who need consistent skin protection and stability during daily use.
- For therapists seeking a versatile cushion that supports postural management and pressure redistribution for clients with complex seating needs.

Key features

- **Firm, pre-contoured base with anatomical well** – promotes stable positioning and effective femoral loading
- **Carveable foam base** – allows for custom build-ups, cut-outs and personalisation
- **Jay Flow fluid tripad** – provides targeted ischial immersion and envelopment for skin protection
- **Compatible with a wide range of accessories** – supports individual postural and positioning needs
- **Air Exchange cover included** – breathable and moisture-wicking, with optional incontinence-resistant cover and solid seat compatibility

ADDITIONAL INFORMATION

Sizes	16" x 16", 16" x 18", 18" x 16", 18" x 18", 18" x 20", 20" x 18"
Brand	Jay