


SKU:

## JAY UNION CUSHION

~~\$899.99~~ **\$899.99**





The Jay Union Cushion blends memory foam and fluid layering to provide effective skin protection and postural support. Its contoured foam base and cover system ensure comfort and hygiene for daily use.







Scan QR code to  
view product online



## VARIATIONS

| Image  | SKU     | Price    | Description | Sizes     |
|--|---------|----------|-------------|-----------|
|    | S14503  | \$659.99 |             | 16" x 18" |
|    | S311285 | \$659.99 |             | 16" x 20" |
|  | S19419  | \$849.99 |             | 18" x 20" |
|  | S13588  | \$849.99 |             | 20" x 18" |

| Image   | SKU    | Price    | Description | Sizes     |
|---|--------|----------|-------------|-----------|
|    | S17284 | \$849.99 |             | 20" x 20" |
|    | S17482 | \$899.99 |             | 22" x 20" |
|  | S12095 | \$729.99 |             | 16" x 16" |
|  | S2124  | \$729.99 |             | 18" x 16" |

| Image  | SKU    | Price    | Description | Sizes     |
|--|--------|----------|-------------|-----------|
|  | S10460 | \$729.99 |             | 18" x 18" |
|  | S2037  | \$729.99 |             | 20" x 16" |

## PRODUCT DESCRIPTION

The Jay Union Cushion is designed for individuals at moderate to high risk of pressure injury who need reliable comfort and positioning during daily activities. Its multi-layered construction includes a gel-infused Visco memory foam top layer that promotes immersion, comfort and heat management. Beneath the memory foam is a Jay Flow fluid insert that conforms to the person's body shape, providing consistent pressure redistribution and reducing the risk of skin breakdown. The contoured foam base features lateral pelvic supports, enhancing stability and optimal alignment. The dual-cover system includes a moisture-resistant inner cover with an Aqua-guard zipper and an antimicrobial outer cover made with silver-impregnated Ionic fibre. This design supports hygiene and reduces odours, making the Jay Union Cushion easy to maintain.

### Who it's for

- For clients at moderate to high risk of pressure injury who need enhanced comfort and postural support during prolonged sitting.
- For therapists seeking a versatile cushion that balances pressure management, stability and hygiene for clients with complex seating needs.

### Key features

- **Gel-infused memory foam top layer** – enhances comfort and reduces heat build-up
- **Jay Flow fluid insert** – conforms to each person's shape for effective pressure redistribution
- **Contoured foam base with lateral supports** – promotes pelvic stability and optimal alignment
- **Moisture-resistant inner cover** – protects foam from fluid ingress and extends cushion lifespan
- **Antimicrobial outer cover with Ionic fibre** – supports hygiene and reduces odours

## ADDITIONAL INFORMATION

|                           |                 |
|---------------------------|-----------------|
| <b>Product Weight Max</b> | 2 kg            |
| <b>Max User Weight</b>    | 136 Kg - 225 Kg |
| <b>Brand</b>              | Jay             |