



SKU:

PRIDE MOBILITY JAZZY AIR 2.0 POWER CHAIR

\$7,999.99

The Pride Mobility Jazzy Air 2.0 Powered Wheelchair features a seat lift that rises 30.4 cm in 11 seconds. It supports indoor and outdoor use with advanced suspension and mid-wheel drive.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Color
	S12452	\$7,999.99		White
	S304594	\$7,999.99		Black

PRODUCT DESCRIPTION

The Pride Mobility Jazzy Air 2.0 Powered Wheelchair transforms the powerchair experience with a rapid-rise seat lift that allows face-to-face social interaction and access to everyday spaces. At the touch of a button, the seat can be raised or lowered in just seconds.

The Jazzy Air's mid-wheel drive and Active-Trac Suspension support a smooth ride on varied terrain, while LED lighting enhances safety.

The compact narrow base and redesigned lift function improves indoor navigation and compatibility with standard height tables or desks. Armrests, foot platform and headrest are all highly adjustable, and under-seat storage offers practical convenience.

Who it's for

- For clients who benefit from powered elevation for interaction and access.
- For therapists prescribing a height-adjustable solution for mobility and participation.

Key features

- **30.4 cm seat elevation in 11 seconds** – encourages social interaction and access
- **Drivable while elevated** – maintain movement at height
- **Active-Trac Suspension & Mid-Wheel 6 Technology** – stability indoors and outdoors
- **Redesigned seat-to-floor height** – enabling user to sit comfortably at standard desks and tables
- **Under-seat storage** – practical for daily essentials
- **Bright, long lasting LED front and rear lights** – enhances visibility and safety
- **Foot platform with 5 cm adjustment** – supports elevated and seated positions
- **Adjustable armrests and headrest** – for personalised comfort

ADDITIONAL INFORMATION

Brand	Pride Mobility
Max SWL	136 kg (300 lbs.)
Maximum Speed	6.4 km/h (4 mph)
Ground Clearance	2.54 cm (1") at motor bracket
High-Back Seat to Floor Range	45.7 x 45.7 – 50.8 cm (18" x 18" – 20")
Turning Radius	61.5 cm (24.25")
Overall Size	Length: 106.7 cm (42") Width: 64.1 cm (25.25")
Drive Wheels	25.4 cm (10") flat free
Caster Wheels	15.24 cm (6") flat free
Suspension	Active-Trac (ATX) with Mid-Wheel 6® Technology
Range	Up to 27 km (16.8 miles)
Component Weights	Base: 40 AH Batteries 120 kg 40 AH Batteries: 14.28 kg
Drivetrain	Two motor, Mid-Wheel 6®
Batteries	(2) 12V, 40 AH