



SKU:


## **K-DEE II HIGH-LOW STANDARD HOSPITAL BED**

**~~\$3,699.99~~ \$3,699.99**

The K-Dee II Standard Bed is the second generation of the popular K-Dee line. These beds are fully featured with five-way adjustment, including height, backrest, knee break, Trendelenburg and reverse Trendelenburg. They are compact and lightweight and easy to assemble and disassemble. This bed has optional accessories including full length or three-quarter side rails for added safety.

Available in both Single and King Single sizes. For taller users an optional length extension can be added to increase the length by 200 mm.

The K-Dee II is designed to be low-maintenance, reliable and robust.



Scan QR code to  
view product online

VARIATIONS

| Image                                                                            | SKU    | Price      | Description | Size        |
|----------------------------------------------------------------------------------|--------|------------|-------------|-------------|
|  | S10209 | \$2,899.99 |             | Single      |
|  | S12391 | \$3,699.99 |             | King Single |

## **PRODUCT DESCRIPTION**

### **Key Features**

- Single or King Single sizes available
- Beech wood grain head and foot board panel
- Can be fully disassembled for transport, and stored in a secured position
- Handset with lockout function
- Trendelenburg and reverse Trendelenburg functionality
- 4 section power adjustable bed
- Height adjustable to suit user needs
- Back rest powered tilt adjustable
- Knee break powered adjustable
- Enclosed individual locking castors
- AkzoNobel quality powder coat
- Dewert actuators
- E1 and class 3 fire retardant
- Tested to: BS EN 60601-1:2006 and BS EN 60601-2-252:2010 Safety Standard

### **Options and Accessories**

- 200 mm length extension option
- Rails: full length steel rails, three-quarter rails

## ADDITIONAL INFORMATION

|                                 |                                                                                                                                                                                                  |
|---------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Maximum user weight</b>      | 213Kg                                                                                                                                                                                            |
| <b>Overall Length</b>           | 2100mm, 2300mm with extension                                                                                                                                                                    |
| <b>Overall Width</b>            | (S) 980mm, (KS) 1160mm                                                                                                                                                                           |
| <b>Mattress platform length</b> | 1900mm, 2100mm with extension                                                                                                                                                                    |
| <b>Mattress platform width</b>  | (S) 870mm, (KS) 1050mm                                                                                                                                                                           |
| <b>Mattress platform height</b> | 335 - 725mm                                                                                                                                                                                      |
| <b>Backrest adjustment</b>      | 0 - 70 °                                                                                                                                                                                         |
| <b>Knee break adjustment</b>    | 0 - 23 °                                                                                                                                                                                         |
| <b>Leg rest adjustment</b>      | -27 - 0 °                                                                                                                                                                                        |
| <b>Size</b>                     | Single, King Single<br>- Bed frame: 5 years<br>- Electrical components (including motors, excluding handset and batteries): 2 years<br>- Castors: 2 years<br>- Handsets and batteries: 12 months |
| <b>Warranty</b>                 |                                                                                                                                                                                                  |