



SKU:


K-DEE II HIGH-LOW STANDARD HOSPITAL BED

~~\$3,699.99~~ **\$3,699.99**

The K-Dee II Standard Bed is the second generation of the popular K-Dee line. These beds are fully featured with five-way adjustment, including height, backrest, knee break, Trendelenburg and reverse Trendelenburg. They are compact and lightweight and easy to assemble and disassemble. This bed has optional accessories including full length or three-quarter side rails for added safety.



Available in both Single and King Single sizes. For taller users an optional length extension can be added to increase the length by 200 mm.

The K-Dee II is designed to be low-maintenance, reliable and robust.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S10209	\$2,899.99		Single
	S12391	\$3,699.99		King Single

PRODUCT DESCRIPTION

Key Features

- Single or King Single sizes available
- Beech wood grain head and foot board panel
- Can be fully disassembled for transport, and stored in a secured position
- Handset with lockout function
- Trendelenburg and reverse Trendelenburg functionality
- 4 section power adjustable bed
- Height adjustable to suit user needs
- Back rest powered tilt adjustable
- Knee break powered adjustable
- Enclosed individual locking castors
- AkzoNobel quality powder coat
- Dewert actuators
- E1 and class 3 fire retardant
- Tested to: BS EN 60601-1:2006 and BS EN 60601-2-252:2010 Safety Standard

Options and Accessories

- 200 mm length extension option
- Rails: full length steel rails, three-quarter rails

ADDITIONAL INFORMATION

Maximum user weight	213Kg
Overall Length	2100mm, 2300mm with extension
Overall Width	(S) 980mm, (KS) 1160mm
Mattress platform length	1900mm, 2100mm with extension
Mattress platform width	(S) 870mm, (KS) 1050mm
Mattress platform height	335 - 725mm
Backrest adjustment	0 - 70 °
Knee break adjustment	0 - 23 °
Leg rest adjustment	-27 - 0 °
Size	Single, King Single - Bed frame: 5 years - Electrical components (including motors, excluding handset and batteries): 2 years - Castors: 2 years - Handsets and batteries: 12 months
Warranty	