




SKU:

KAYE THERAPY BENCH

The Kaye Therapy Bench is ideal for home, therapy, or school programs and for children who do not need back support in sitting but need a firm seat that can be adjusted to ensure the child's feet rest on the floor. Manufactured by Kaye Products, the bench can be used for learning, playing and daily living.



Scan QR code to
view product online

PRODUCT DESCRIPTION

Benches can be used in a variety of therapy settings, including physical, occupational and speech therapy, particularly when the child needs a stable sitting surface to achieve maximum arm and hand movements.

The solid, stable construction makes the Kaye Therapy Bench ideal for children learning to move from one level to another, such as moving from the floor to sitting, from sitting to standing, or from a seated position to their walker.

It is quick and easy to raise, lower or tilt the seat surface due to the slotted legs. Many children benefit from a tilted surface, which places more weight on their feet and aids postural control for better arm and hand function.

- An adjustable bench ideal for home, therapy, or school programs
- Available in five sizes

ADDITIONAL INFORMATION

Max weight

105 kg