

SKU: S16271



KENNINGTON TWIN SEAT DUAL MOTOR LIFT CHAIR WITH HEADREST AND LUMBAR - HOME THEATER RECLINERS

\$5,999.99

Sit back, relax and relish Kennington's top of the line duo Lift Chair. The recliner offers supreme comfort, a touch of class and an easy electric recline. This one of a kind dual motor recliner features an independent head and leg rest function, cup holders, handwand with a USB port and console storage.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Kennington Twin Seat Dual Motor lift chair provides great relaxation, a hint of elegance, and a simple motorized recline. This one-of-a-kind dual motor recliner includes cup holders, a handwand with a USB port, and console storage in addition to an independent head and leg rest function.

Single vs twin vs quad motor lift chairs - what's the difference?

All ILS lift chairs have a lifting action, to support the user to stand from a seated position.

- Single motor chairs move the back rest and leg rest at the same time.
- Twin motor chairs allow the backrest and leg rest to move independently of each other, typically allowing for a greater degree of leg elevation and greater degree of recline compared to a standard single motor chair.
- Quad motor chairs allow independent movement of headrest, lumbar, backrest and footrest allowing for the greatest level of re-positioning of all chairs.

ADDITIONAL INFORMATION

| | |
|----------------------|----------------------------|
| Max SWL | Each Chair 158.7kg Max SWL |
| Seat Width | 57.1cm |
| Seat Depth | 54.6cm |
| Seat-to-Floor | 49.5cm |
| Chair Width | 193cm |
| Chair Depth | 106.7cm |
| Chair Height | 106.7cm |
| Colour | Stonewash (Gun Metal) |
| Brand | Theorem |
| Warranty | 2 years |