



**SKU:**

## **KINGSTON SWIVEL LOWBACK DAY CHAIR**

**\$899.99**

The Kingston Swivel Low Back Day Chair is a versatile mobile seating option designed to assist with dining transfers and positioning. It combines adjustable support with smooth mobility to suit various care environments.



Scan QR code to  
view product online

## VARIATIONS

Image	SKU	Price	Description	Style
	S10016	\$899.99		Beige Vinyl
	S10074	\$899.99		Grey Light Fabric

## PRODUCT DESCRIPTION

The Kingston Swivel Low Back Day Chair is a supportive and manoeuvrable seating option for those requiring mobility assistance in dining and communal settings.

This low back chair is fitted with four castors – two of which are lockable – allowing safe, controlled movement. A rear push handle enhances caregiver control during repositioning.

The chair includes height-adjustable legs and drop-down armrests to accommodate individual needs and make side transfers more manageable.

Available in both fabric and vinyl upholstery, the Kingston offers easy-care surfaces suited to clinical or home environments.

Its padded seat, backrest and arm pads provide additional comfort for extended use.

### Who it's for

- For clients who need mobile seating with adjustable positioning for easier table access.
- For therapists seeking a functional, lightweight day chair to support seated transfers.

### Key features

- **Height-adjustable legs** – support correct seating posture and comfort
- **Drop-down armrests** – assist with safe and easy side transfers
- **Four castors (two lockable)** – enable smooth repositioning and stability when needed
- **Rear push handle** – improves manoeuvrability and carer control
- **Padded seat, backrest and arms** – enhance long-term sitting comfort
- **Vinyl or fabric upholstery options** – support hygienic maintenance and personal preference

## **ADDITIONAL INFORMATION**

<b>Seat Width</b>	46 cm
<b>Seat Depth</b>	46 cm
<b>Floor To Seat Height Adjustment</b>	45 cm – 57.5 cm
<b>Arm Rest Adjustable Height</b>	0 cm – 26 cm
<b>Max Weight Capacity / SWL</b>	160 kg