



**SKU:** S598

## **HEMECRAFT KITCHEN PERCHING STOOL**

**\$285.99**

The Homecraft Kitchen Perching Stool provides seated support at bench height, helping reduce strain on the back and legs. Its adjustable legs padded seat and stable frame promote safety and comfort during use in the kitchen or household.



Scan QR code to  
view product online

## PRODUCT DESCRIPTION

The Homecraft Kitchen Perching Stool is designed to assist people who tire easily when standing for daily tasks. It features a broad, forward-sloping padded seat and supportive backrest, both upholstered in wipe-clean vinyl for hygiene and comfort.

Height-adjustable legs allow the stool to be set at a safe and practical height for each individual. The lightweight aluminium frame is rust-resistant and easy to reposition, while large non-slip rubber feet offer a stable base across a range of floor surfaces.

Ideal for kitchen tasks or use in other standing environments, the stool helps improve posture and reduce fatigue for those with limited strength or mobility.

### Who it's for

- For clients who require seated support to reduce fatigue while preparing meals or managing household tasks.
- For therapists supporting people with mobility limitations to remain independent during daily activities.

### Key features

- **Height-adjustable legs** – accommodates different individuals and benchtop heights
- **Padded vinyl seat and backrest** – sloped seat supports posture and reduces fatigue
- **Lightweight aluminium frame** – easy to move and rust-resistant for long-term use
- **Large non-slip feet** – enhance safety on kitchen and household surfaces
- **Durable construction** – supports hygiene, comfort and stability for regular use

## **ADDITIONAL INFORMATION**

<b>Maximum User Weight</b>	120 kg
<b>Width Between Arms</b>	45.5 cm
<b>Seat Depth</b>	30 cm
<b>Seat To Top of Backrest</b>	30 cm
<b>Leg Height Adjustment</b>	50 cm
<b>Floor To Seat</b>	63.5 cm