


**SKU:** S12286

## **HERO MEDICAL KNEE WALKER - LIGHTWEIGHT**

**\$349.99**

The lightweight Hero Medical knee walker only weighs 9.1kg and can fold away for easy storage. Hero Knee walker has a durable steel frame with a safe working load of 136kg with 8" PVC wheels so you won't need to worry about it going flat, and a comfortable polyurethane knee pad. A handlebar brake will allow the user to be in control for maximum security. Buy with confidence knowing that the product has a 12-month warranty.



Scan QR code to  
view product online



## **PRODUCT DESCRIPTION**

The Hero medical Knee Walker has a durable steel frame with a safe working load of 136KG, 8" PVC wheels that won't go flat, and a comfortable polyurethane knee pad. A handlebar brake puts the user in control for maximum security.

## ADDITIONAL INFORMATION

|                         |              |
|-------------------------|--------------|
| <b>Brand</b>            | Hero Medical |
| <b>Max SWL</b>          | 136kg        |
| <b>Weight</b>           | 9.1kg        |
| <b>Handlebar height</b> | 860-1010mm   |
| <b>Knee pad height</b>  | 490mm-590mm  |
| <b>Frame</b>            | Steel        |
| <b>Wheels</b>           | PVC          |
| <b>Knee pad</b>         | Polyurethane |