




SKU:

## **DAYS HEALTHCARE LINK WHEELCHAIR**

**~~\$1,299.99~~ \$1,299.99**

The Days Healthcare Link Wheelchair is a lightweight adaptable unit and is suitable for full-time users and is designed for maximum practicality and adjustability. For comfort, safety and convenience, the Link is fitted with swing-back and detachable armrests, an angle adjustable footplate, and a heel strap. The Link is an attractive and adaptable chair which offers six different seating widths to suit your individual needs.



Scan QR code to  
view product online

VARIATIONS

Image	SKU	Price	Description	Size
	S2512	\$899.99		16 inches
	S3064	\$1,299.99		20 inches
	S2513	\$1,299.99		22 inches

## **PRODUCT DESCRIPTION**

Optional additions such as height-adjustable armrests and elevating leg rests can be added for added comfort. The Days Healthcare Link Wheelchair is equipped with puncture-resistant tyres for further peace of mind while out and about. For added safety on hills, rear anti-tippers are included as standard equipment.

### **Features**

- Quick release wheels
- Solid tyres
- 16", 18", 20", and 22" seat widths available
- Swing in/out leg rests
- Swing back removable armrests
- Height adjustable back upholstery
- Adjustable castors
- Multi position rear wheels
- Drum brakes on the 20"-22" seat widths
- Anti-tipper is included as standard.

## ADDITIONAL INFORMATION

<b>Brand</b>	Days Healthcare
<b>Seat depth</b>	43cm (16.9 inches)
<b>Seat width</b>	16" 40cm (15.7 inches)
	20" 50.8cm (20 inches)
	22" 55.8cm (21.9 inches)
<b>Overall width</b>	62cm (24.4 inches)
<b>Overall height</b>	98cm (38.5 inches)
<b>Max SWL</b>	User Weight 16" & 18" 130kg (286 lbs)
	User Weight 20" & 22" 150kg (330 lbs)
<b>Weight of chair</b>	17.5 kg
<b>Weight without wheels</b>	14kg Warranty 24 months on frame & crossbars. Does not include tyres, upholstery or arm pads.
<b>wheelchair-type</b>	Self Propelled