




SKU:

## **DAYS HEALTHCARE LINK WHEELCHAIR**

**~~\$1,299.99~~ \$1,299.99**

The Days Healthcare Link Wheelchair is a lightweight adaptable unit and is suitable for full-time users and is designed for maximum practicality and adjustability. For comfort, safety and convenience, the Link is fitted with swing-back and detachable armrests, an angle adjustable footplate, and a heel strap. The Link is an attractive and adaptable chair which offers six different seating widths to suit your individual needs.



Scan QR code to  
view product online

VARIATIONS

| Image                                                                            | SKU   | Price      | Description | Size      |
|----------------------------------------------------------------------------------|-------|------------|-------------|-----------|
|  | S2512 | \$899.99   |             | 16 inches |
|  | S3064 | \$1,299.99 |             | 20 inches |
|  | S2513 | \$1,299.99 |             | 22 inches |

## **PRODUCT DESCRIPTION**

Optional additions such as height-adjustable armrests and elevating leg rests can be added for added comfort. The Days Healthcare Link Wheelchair is equipped with puncture-resistant tyres for further peace of mind while out and about. For added safety on hills, rear anti-tippers are included as standard equipment.

### **Features**

- Quick release wheels
- Solid tyres
- 16", 18", 20", and 22" seat widths available
- Swing in/out leg rests
- Swing back removable armrests
- Height adjustable back upholstery
- Adjustable castors
- Multi position rear wheels
- Drum brakes on the 20"-22" seat widths
- Anti-tipper is included as standard.

## ADDITIONAL INFORMATION

|                              |                                                                                               |
|------------------------------|-----------------------------------------------------------------------------------------------|
| <b>Brand</b>                 | Days Healthcare                                                                               |
| <b>Seat depth</b>            | 43cm (16.9 inches)                                                                            |
| <b>Seat width</b>            | 16" 40cm (15.7 inches)                                                                        |
|                              | 20" 50.8cm (20 inches)                                                                        |
|                              | 22" 55.8cm (21.9 inches)                                                                      |
| <b>Overall width</b>         | 62cm (24.4 inches)                                                                            |
| <b>Overall height</b>        | 98cm (38.5 inches)                                                                            |
| <b>Max SWL</b>               | User Weight 16" & 18" 130kg (286 lbs)                                                         |
|                              | User Weight 20" & 22" 150kg (330 lbs)                                                         |
| <b>Weight of chair</b>       | 17.5 kg                                                                                       |
| <b>Weight without wheels</b> | 14kg Warranty 24 months on frame & crossbars. Does not include tyres, upholstery or arm pads. |
| <b>wheelchair-type</b>       | Self Propelled                                                                                |