


SKU: S1658

LONG HANDI EASIREACHER

~~\$65.99~~ **\$65.99**



The Homecraft Pick-Up Reacher is a lightweight, robust reacher ideal for retrieving or reaching various items in the home, garden or stores. It has twin gripping jaws with sculptured liners and a trigger designed for comfort. A magnetic head enables easy retrieval of pins, paper clips and other small, metal items. The handle is large and roomy to allow for the use of all fingers, and the short trigger distance is ideal for those with stiff, arthritic fingers or limited mobility. It has an easily rotatable head and a removable stick clip to connect the reacher to wheelchairs, walking frames and other mobility devices.



Scan QR code to
view product online

VARIATIONS

Image	SKU	Price	Description	Length
	S1043	\$56.99		60 cm
	S1042	\$59.99		75 cm
	S1044	\$65.99		90 cm

PRODUCT DESCRIPTION

The Homecraft Pick-Up Reacher is easy to operate for anyone, especially users with hand and wrist injuries. The rotating head keeps users from twisting their wrists in painful ways, while the handle has four finger grooves, instead of the standard two, for extra grip force.

Key Features

- The trigger is designed for comfort and is large and roomy to allow use of all fingers. The short trigger distance is ideal for stiff, arthritic fingers. The trigger projects through the top of the handle to allow it to be squeezed in the conventional manner or pushed using the palm of the hand.
- The jaws have sculptured liners which give an excellent grip. The Pick-Up is operated by a wire not a cord, eliminating fraying or stretching.
- Magnet attached to the jaw which allows easy retrieval of pins, paper clips and other small ferrous items
- The head can easily rotate to suit individual tasks. The Pick-Up is supplied with a removable stick clip, allowing it to be clipped on to wheelchairs and walking frames etc.

ADDITIONAL INFORMATION

Brand	Performance Health
--------------	--------------------