




**SKU:** S1658

## **LONG HANDI EASIREACHER**




**~~\$65.99~~ \$65.99**

The easireacher adds inches to arm length, so picking up objects doesn't have to include stretching and bending. An economical yet durable reacher constructed with wire-driven jaws, eliminating a string that can stretch and break.



Scan QR code to  
view product online

### VARIATIONS

Image	SKU	Price	Description	Length
	S1043	\$56.99		60 cm
	s1042	\$59.99		75 cm
	S1044	\$65.99		90 cm

## **PRODUCT DESCRIPTION**

The trigger, operated by four fingers instead of the usual two is easy to grip and operate and slip-resistant jaws can pick up even extremely small items. Includes a magnet for picking up metal objects and a hook that aids with dressing. The lightweight frame resists twisting and bending when picking up heavier objects.

## **ADDITIONAL INFORMATION**

**Brand**

ILS ADL