




SKU:

MAGIC MOBILITY EXTREME X8

The Magic Mobility Extreme X8 is an all-terrain powerchair that lets you explore off-road trails with ease.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The Extreme X8 is designed to take on any terrain, from sand dunes to snowy hillsides. Its 4-pole motors deliver high speed and torque, while the low-pressure off-road knobby tyres provide a smooth ride over bumps and uneven surfaces. With a top speed of 10 km/h, the X8 gives you the freedom to explore the great outdoors at your own pace. Proudly made in Australia by Magic Mobility, the Extreme X8 is built to last and is perfect for anyone who wants to venture off the beaten path.

ADDITIONAL INFORMATION

Brand	Magic Mobility
Type	Off Road
Drive Wheel Configuration	All Wheel
Weight Capacity	182 kg, 155 kg (Tilt and/or Lift)
Transit Approved	Yes
Max Speed	10 km/h
Base Width	28"
Base Length	41"
Pediatric Seating	No
Seat Width Range	12" - 24"
Seat Depth Range	12" - 24"
Seat To Floor Height Range	18.5" - 19.5"
Posterior Tilt Option	0° - 45°
Anterior Tilt Option	No
Recline Option	Yes
Legrest Elevation Option	Yes
Seat Elevation Option	12"
Power Standing Option	No