




**SKU:**

## **ROYALE MAYFAIR LUXURY LEATHER LIFT CHAIR (TWIN MOTOR)**

**\$0.00**

Exclusive to ILS, the Royale Mayfair Luxury leather lift chair combines premium materials with quality twin motor functionality, offering individuals a comfortable, supportive and attractive lift chair solution. The twin motor system offers users independent control of the backrest and legrest to ensure flexible and comfortable positioning.



Scan QR code to  
view product online

## VARIATIONS

Image	SKU	Price	Description	Colours
	S13431	\$2,799.99		Charcoal Black
	S13433	\$2,699.99		Taupe Beige
	S13432	\$2,799.99		Espresso Dark Brown

## PRODUCT DESCRIPTION

The Royale Mayfair Luxury Leather Lift Chair (Twin Motor) offers a twin motor design that allows separate adjustment of the backrest and legrest, which provides greater flexibility in seating positions. The lift mechanism quietly assists transitions into standing and sitting positions. Designed with a three-section waterfall backrest and upholstered in premium leather from high-quality hides, this chair delivers breathable comfort and visual appeal.

It's suited to individuals who want both supportive function and contemporary styling.

### Who it's for

- For clients who value premium materials and want greater positioning control.
- For therapists needing a stylish twin motor chair with leather finish to support their client.

### Key features

- **Twin motor functionality** - independent backrest and legrest control
- **Luxury leather upholstery** - durable, breathable and soft to the touch
- **Three-section waterfall backrest** - promotes ergonomic alignment and support
- **Premium foam seat cushioning** - for long-term comfort and postural support
- **Smooth vertical lift** - assists with easy and safe transitions
- **Battery back-up** - ensures reliability during power interruptions
- **Easy-to-use hand control** - intuitive design
- **Side pockets** - practical and accessible storage

## ADDITIONAL INFORMATION

<b>Safe Working Load (SWL)</b>	150 kg
<b>Overall Width</b>	85 cm
<b>Overall Depth</b>	85 cm
<b>Overall Height</b>	105 cm
<b>Seat Width</b>	52 cm
<b>Seat Depth</b>	50 cm
<b>Seat to Floor Height</b>	48 cm
<b>Number of Motors</b>	Twin
<b>Brand</b>	Royale Medical