


SKU: S2646

ANALOG TENS MACHINE BY ILS PHYSIO

~~\$89.99~~ \$62.99



A regular TENS with adjustable pulse width, pulse rate, intensity, modes and an accurate timer of 15, 30 and 60 minutes, will stop output at preset time. The unit comes with the TENS stimulator unit, carrying case, 2X TENS leads, 9V battery, operating manual and packet of 4 Electrodes. Includes a belt clip.



Scan QR code to
view product online

PRODUCT DESCRIPTION

The ProTENS unit has a fully adjustable parameter and provides relief from acute or chronic pain such as back pain, sports injury, general muscle pain, etc. Battery powered electrical unit which uses electrodes to deliver electrical pulses to the nerve fibres relieving pain.

ADDITIONAL INFORMATION

| | |
|------------------------|---|
| Brand | ILS physio |
| Channel | Dual, isolated between channels |
| Modes | Continuous, Burst, Modulation |
| Pulse Intensity | Adjustable 0-80mA peak into 500 ohm load each channel, constant current |
| Pulse Rate | 2Hz-150Hz (adjustable) |
| Pulse Width | 30µs-260µs (adjustable) |
| Time | Continuous, 30 min., 60 min. |
| Burst Mode | Burst consists 2 burst per sec at 100 Hz |
| Wave Form | Asymmetrical Bi-Phasic square pulse |
| Dimensions | 95(H) x 60(W) x 23(D) mm |
| Weight | 115 grams (battery included) |
| Warranty | 3 Years |